

# Data Review and Analysis with iC IR™



iC IR displays dynamic linked views in live and review experiments. Four main 'views' provide a cross-section of experiment information, and a selection or change in one view dynamically appears in the other linked views. A 'toolbox' **I** enables you to control a variety of tasks.

## I. Linked Views

- **Trends Viewer** **A**  
Displays real-time trends of each user-defined peak profile and ConcIRT™ LIVE! component in the experiment. Trends can also include Quant and Referee components plus user-defined math calculations.
- **Data/Events Viewer** **B**  
Select from two tabs. Data tab shows the time and value of each sample collected. Events tab shows a record of all user annotations and system/audit messages.
- **Spectra Viewer** **C**  
Shows a 2D view of pinned and selected spectra and allows definition of peak profiles.
- **Surface Viewer** **D**  
Presents a 3D view of the entire experiment and gives you the following options:
  - **Rotate or move**—Place cursor on the surface view, click the left mouse button and drag to ideal rotation. Use the right mouse button to move the 3D view.
  - **Enlarge or reduce size**—Place cursor in the view and rotate the mouse wheel.
  - **Additional options**—Right-click anywhere in the view to see more options.

## II. Live Experiment Toolbar **E**

- Check current status of experiment by noting the background color of the toolbar:  
Green=Running, Yellow=Paused, Blue=Completed, Red=Error condition

- **F** Control experiment status using the **Play/Pause** and **Stop** buttons.  
▶ = Play    ⏸ = Pause    ⏹ = Stop
- **G** Click the **Schedule** button to add/ edit/delete a phase of the experiment.
- **K** Click the **Add Referee** button after taking a manual sample for HPLC or another chromatography method.

## III. Annotations

Make any number of annotations on conditions or actions such as the beginning or end of additions, when a sample is extracted for offline analysis, modifications in reaction parameters, reaction color, bubbles, and so on. This provides useful details when you analyze results.

### Adding an Annotation—Three Options

- **H In Live Experiment**—Type the annotation in the Live Experiment Toolbar while the experiment is running. The Events Viewer logs the note.  
To see annotations in the Trend Viewer, click the **Show/Hide Annotation** button .  
To add annotations, click the **Add Annotation** button .
- **A In Trends Viewer**  
Click to set the point where you want to add the annotation and use either the toolbar button or the right-click menu to open the annotation text box.
- **B In Events Viewer**  
Select the Events tab, right-click a row in the events list, and select **Annotation/Add**.  
To edit an annotation, double-click the annotation icon in the Trends or Events viewers or the row in the Events Viewer.

## IV. Peak Profiles and Labels

- **C** In the **Spectra Viewer**, you can define and edit profiles for peaks of interest and create or edit peak labels.

### Defining a Peak Profile—Two Options

- Both options create a peak profile definition that follows the maximum absorbance value over the specified region, based on the currently selected spectrum.
- Double-click above the spectrum to define a typically narrower region centered on the maximum absorbance value of the component of interest.
  - Double-click below the spectrum to define a typically broader region determined by the inflection points of the component of interest.


### Changing a Peak Profile

- Ensure the Show Details option is checked in the context menu (Select **Show Details** button .
- Click the **Peaks** tab along the left side of the details table.
- Double-click a specific Peak Profile or graphically edit a peak by selecting it; then moving it or adjusting the width.
- Select the peak category from the drop-down menu in the Group field.
- Select the peak method from the drop-down menu in the Type field, and click **Apply**.



## Creating/Editing a Peak Label

Create and edit peak labels in the Spectra View as well as in a Spectra Library.

1. Pin a spectrum by selecting it and clicking the **Pin Sample** button.
2. Click the **Add Peak Labels** toolbar button.



3. Click and drag across one or more peaks to select the ones to be labeled. The system creates a label with the peak center point value. Click and drag the label to move it.
4. Double-click on the label to edit it.



To edit an existing label, make sure the **Show Point Labels** button is selected, then double-click the label.

## V. Process Spectral Data

From the Toolbox, select the **Data Treatments** task pane **I**. In the default view, the most commonly used options appear under the following categories:

### Reference Spectra

- **Correct for water vapor**—Apply a correction for water vapor based on a sample collection taken before the experiment (optional for sealed instruments: ReactIR 15, 247, 45P).
- **Subtract spectrum**—Subtract the selected reference spectrum from all sample spectra in the experiment.
- **Load or Remove**—Add reference spectra to or delete from the experiment.

### Functions

The more common functions include:

- **Baseline offset**—Vertically shift each spectrum in the data set to zero at the selected reference wavenumber.
- **Normalize**—Correct each spectrum in the data set to a range of 0 at a single baseline point to 1.0 at the selected reference wavenumber. This normalization is useful to isolate effects such as temperature when you have a known standard in your reaction that should result in a constant peak height.
- **2nd derivative**—Apply a 2nd derivative function to each spectrum in a data set. This produces a spectrum that more easily processes weak remote or overlapping IR bands. The function multiplies results by  $-1$  so traditionally negative peaks appear positive.

### ConcIRT

- **Enable (ConcIRT LIVE!)**—Extract relative concentration profiles for products, intermediates, and starting materials, especially for chemical species that have overlapping peaks.
- **Show Components**—Show component spectra identified by ConcIRT LIVE!

### Display Options

Use the 'Hide the DiComp "Diamond Region" check box to suppress display of the spectral region from 1950 to 2250 $\text{cm}^{-1}$  where the DiComp sensor absorbs the infrared energy. To specify a region other than the default, click the options button.

### Background Replacement

To replace the background with one loaded from an external source, mark "Show all options."

## VI. Spectra Libraries

Create standard chemical references separately from an experiment and store them in spectra libraries.


1. Click **New Spectra Library** from the Start Page.
2. Click either the **Sample Using Wizard**, or **Quick Collect** button to collect samples. Follow the steps under "Reference Spectra Needed?" in the 'Experiment Setup in iCIR' quick reference document.

### Moving Spectra To or From a Library

Add selected spectra to the current experiment from a library, or add selected spectra from an open library to an experiment. Select spectra and use the drag-and-drop method.

## VII. Result Sets

A result set is a separate file that contains trends. Data from a series of related experiments can be placed in one result set for comparison. To create the file, click **New Result Set** from the Start Page, select the trends, and save them using the 'File/Save as...' command.



Add trends to a result set by either of the following methods:

- **J** **In the graph area of the Results Set window**—Right-click on the graph and select **Add Trend** option.
- **A** **In the Trends Viewer for the experiment**—Click to highlight a trend., then drag-and-drop it to the target Result Set tab.

## VIII. Templates

To enable running the same reaction multiple times with different conditions, create templates of frequently used trends and data treatments to reuse in subsequent experiments.

**Watch the "Guided Tour" from the iCIR Start Page to see selected data review and analysis features in action.**

**Mettler-Toledo AutoChem, Inc.**

7075 Samuel Morse Drive  
Columbia, MD 21046 USA  
Telephone +1 410 910 8500  
Subject to technical changes  
P/N: 14690102  
© 06/2011 MK-PB-0012-AC Rev C DCN 559

**METTLER TOLEDO**

[www.mt.com/icir](http://www.mt.com/icir)