

# DIFFRAC.SUITE

- Tutorial

DIFFRAC.EVALUATION PACKAGE  
DIFFRAC.EVA

Original Instructions

The reproduction, transmission or use of this document or its contents is not permitted without express written authority. Offenders will be liable for damages. All rights reserved.

We have checked the contents of this manual for agreement with the hardware and software described. Since deviations cannot be precluded entirely, we cannot guarantee full agreement. However, the data in this manual are reviewed regularly and any necessary corrections are included in subsequent editions. Suggestions for improvement are welcome.

All configurations and specifications are subject to change without notice.

Order no. DOC-M88-EXX201 V4. Updated: Nov 27, 2014.

© 2010 - 2014 Bruker AXS GmbH, Karlsruhe, Germany.

All trademarks and registered trademarks are the sole property of their respective owners.

Printed in the Federal Republic of Germany.

Bruker AXS GmbH  
Östliche Rheinbrückenstr. 49  
76187 Karlsruhe, Germany  
Tel. +49 721 50997-0  
Fax +49 721 50997-5654  
info@bruker-axs.de  
www.bruker.com

# DIFFRAC.EVA Tutorial

## Table of Contents

<b>EVA Tutorial.....</b>	<b>1</b>
<b>Performing a Search/Match Operation.....</b>	<b>1</b>
Step 1: Creating a New EVA Document and Importing Bchips.RAW .....	2
Step 2: Setting the Search Parameters .....	3
Step 3: Matching Procedure .....	5
Step 4: Preparing the Residual Scan .....	6
Step 5: Performing a Subsequent Search on the Residual Scan .....	7
Step 6: Matching Procedure .....	8
<b>Performing an Automatic Search/Match Operation.....</b>	<b>10</b>
Case #1: m1.RAW.....	11
Step 1: Creating a New EVA Document and Importing m1.RAW .....	11
Step 2: Setting the Search Parameters .....	12
Step 3: Automatic Matching Procedure .....	15
Case #2: CPD-2.RAW.....	18
Step 1: Creating a New EVA Document and Importing CPD-2.RAW .....	18
Step 2: Setting the Search Parameters .....	19
Step 3: Automatic Matching Procedure .....	22
Case #3: Mineral mixture .....	23
Step 1: Creating a New EVA Document and Importing Mineral_mixture.brml.....	23
Step 2: Setting the Search Parameters .....	24
Step 3: Automatic Matching Procedure .....	27
Step 4: Saving .....	28
<b>Creating and Using Filter Lists .....</b>	<b>29</b>
Step 1: Creating a New EVA Document and Importing BX100.RAW .....	30
Step 2: Creating a Filter List.....	31
Step 3: Using a Filter List .....	38
<b>Creating a User Database .....</b>	<b>41</b>
Step 1: Creating a New EVA Document and Importing m1.RAW.....	42
Step 2: Performing the Automatic Matching Procedure .....	43
Step 3: Performing the Tune Cell Operation .....	47
Step 4: Creating the User Database.....	48
Step 5: Adding the Calcite pattern to the database.....	50
Step 6: Creating a New Document and Importing LaB6.BRML .....	51
Step 7: Adding the LaB6 DIF to the User Database.....	52
Step 8: Using the User Database during a Search/Match operation .....	55
Step 9: Saving .....	57
<b>Performing a Peak Search.....</b>	<b>58</b>
Step 1: Creating a New EVA Document and Importing Csand.RAW.....	59
Step 2: Setting the Peak Search Parameters.....	60
Step 3: Appending the Peaks .....	61
Step 4: Saving .....	62
<b>Computing <math>K\alpha_2</math> Stripping .....</b>	<b>63</b>
Step 1: Creating a New EVA Document and Importing Quartz.RAW .....	64
Step 2: Computing the $K\alpha_2$ Stripping.....	65
Step 3: Saving .....	66
<b>Smoothing Scans .....</b>	<b>67</b>
Step 1: Creating a New EVA Document and Importing Quartz.RAW .....	68

Step 2: Smoothing the Scan .....	69
Step 3: Saving .....	72
<b>Removing Aberrant Points .....</b>	<b>73</b>
Step 1: Creating a New EVA Document and Importing Spikes.RAW .....	74
Step 2: Removing Aberrant Points .....	75
Step 3: Saving .....	76
<b>Computing Areas .....</b>	<b>77</b>
Step 1: Creating a New EVA Document and Importing ET20.RAW .....	78
Step 2: Selecting and Computing an Area .....	79
Step 3: Saving .....	81
<b>Adding Scans .....</b>	<b>82</b>
Step 1: Creating a New EVA Document and Importing ET20.RAW, ET50.RAW and ET80.RAW .....	83
Step 2: Adding the Scans .....	84
Step 3: Saving .....	85
<b>Subtracting Scans .....</b>	<b>86</b>
Step 1: Creating a New EVA Document and Importing ET20.RAW and ET50.RAW .....	87
Step 2: Subtracting the Scans .....	88
Step 3: Saving .....	89
<b>Merging Scans .....</b>	<b>90</b>
Step 1: Creating a New EVA Document and Importing Corundum-1.RAW, Corundum-2.RAW and Corundum-3.RAW .....	91
Step 2: Merging the Scans .....	92
Step 3: Saving .....	93
<b>Normalizing Scans .....</b>	<b>94</b>
Step 1: Creating a New EVA Document and Importing Znox1.RAW and Znox2.RAW .....	95
Step 2: Normalizing the Scans .....	96
Step 3: Saving .....	97
<b>Computing the Crystallinity .....</b>	<b>98</b>
Step 1: Creating a New EVA Document and Importing m2.RAW .....	99
Step 2: Computing the Crystallinity .....	100
Step 3: Checking the Results .....	101
Step 4: Saving .....	102
<b>Simulating a Slit Mode .....</b>	<b>103</b>
Step 1: Creating a New EVA Document and Importing LaB6.BRML .....	104
Step 2: Performing the Search/Match Operation .....	105
Step 3: Simulating a Slit Mode .....	106
Step 4: Saving .....	107
<b>Performing the Semi-Quantitative Phase-Analysis and Comparing with a Chemical Analysis .....</b>	<b>108</b>
Step 1: Creating a New EVA Document and Importing BX100.RAW .....	109
Step 2: Performing the Search/Match Operation .....	110
Step 3: Performing the Semi-Quantitative Phase Analysis .....	113
Step 4: Saving .....	115
<b>Using the d Multiplied By Tool.....</b>	<b>116</b>
Step 1: Creating a New EVA Document and Importing Csand.RAW .....	117
Step 2: Performing the Search/Match Operation .....	118
Step 3: Using the d Multiplied By Tool.....	120
Step 4: Saving .....	121
<b>Using the Tune Cell Tool .....</b>	<b>122</b>
Step 1: Creating a New EVA Document and Importing Francolite.RAW .....	123
Step 2: Performing the Search/Match Operation .....	124
Step 3: Performing the Tune Cell Operation .....	125
Step 4: Saving .....	126
<b>Displaying a Pattern with Another Wavelength .....</b>	<b>127</b>
Step 1: Creating a New EVA Document and Importing Francolite.RAW .....	128
Step 2: Performing the Search/Match Operation .....	129
Step 3: Changing the Wavelength.....	130
Step 4: Saving .....	131
<b>Working with the Waterfall View .....</b>	<b>132</b>

Step 1:	Creating a New EVA Document and Importing Ep333f.RAW.....	133
Step 2:	Creating and Working with the Waterfall View.....	134
Step 3:	Saving .....	136
<b>Working with the 2D View .....</b>		<b>137</b>
Step 1:	Creating a New EVA Document and Importing Ep333f.RAW.....	138
Step 2:	Creating and working with the 2D view.....	139
Step 3:	Saving .....	141
<b>Sorting a Multi-range Scan by Temperature in a 2D View .....</b>		<b>142</b>
Step 1:	Creating a New EVA Document and Importing Guil4.RAW.....	143
Step 2:	Creating the 2D View with Levels .....	144
Step 3:	Sorting the Scans on the Y-Axis by Temperature.....	145
Step 4:	Saving .....	146
<b>Working with the PIP and VIP Views .....</b>		<b>147</b>
Step 1:	Creating a New EVA Document and Importing m1.RAW .....	148
Step 2:	Creating and working with the PIP and VIP views.....	149
Step 3:	Saving .....	152
<b>Creating a Label from a Peak.....</b>		<b>153</b>
<b>Loading and Integrating Mergeable 2D Frames.....</b>		<b>156</b>
Step 1:	Creating a New EVA Document and Importing the 2D Frames.....	157
Step 2:	Selecting the Integration Cursor and Drawing the Integration Area .....	158
Step 3:	Displaying the Cursors Preview Tool and Integration.....	159
<b>Loading and Integrating Stackable 2D Frames.....</b>		<b>160</b>
Step 1:	Creating a New EVA Document and Importing the 2D Frames.....	161
Step 2:	Selecting the Integration Cursor and Drawing the Integration Area .....	162
Step 3:	Displaying the Cursors Preview Tool and Integrate.....	163
<b>Rocking Curve Integration on Stackable 2D Frames .....</b>		<b>164</b>
Step 1:	Creating a New EVA Document and Importing the 2D Frames.....	165
Step 2:	Selecting the Integration Cursor and Drawing the Integration Area .....	166
Step 3:	Displaying the Cursors Preview Tool and Integrate.....	167
<b>Adjusting the Default Mask of a 2D Frame .....</b>		<b>168</b>
Step 1:	Creating a New EVA Document and Importing the 2D Frames.....	169
Step 2:	Check the Integration with the Full Frame Cursor .....	170
Step 3:	Changing the Default Mask to Exclude Non-Exposed Areas .....	171
Step 4:	Selecting the Integration Cursor and Integration .....	172
<b>Using Printing Options .....</b>		<b>173</b>
	Vertical Table .....	173
	Repeat Columns.....	175

(intentionally left blank)

## EVA Tutorial

### Performing a Search/Match Operation



#### NOTE

A reference database is required to perform a Search/Match operation.

This tutorial was prepared using PDF 4+ 2009 as a reference database. You may have another reference database.

If you have no database, you cannot perform this part of the tutorial concerning the Search/Match.

The following seven-step procedure describes the basic procedure for running EVA Search/Match, using the default parameters.

The scan used is held as a tutorial file, BChips.RAW, found in the Tutorial directory.

These are boiling chips coming from the stockroom of the Department of Chemistry, North Dakota State University.

The boiling chips are crushed Dolomite/calcite marble. Quartz from the precursor limestone is a third, routinely identified phase. In the metamorphic rock, weak peaks from one or more layer Silicates are detected. They may be too weak to identify by computer routines, but important low angle *d*-spacings, familiar to clay mineralogists, suggest a Chlorite (e.g. Clinochlore) phase.

#### Steps

1. Creating a new EVA document and importing BChips.RAW.
2. Setting of the search parameters and performing the initial search.
3. Matching procedure:  
Identification of Dolomite (pattern 00-036-0426), Calcite (01-89-1304) and Quartz (01-075-8322)
4. Preparing the residual scan.
5. Performing a subsequent search on the residual scan.
6. Matching procedure  
Identification of Clinochlore (01-083-1365)
7. Saving the EVA document containing the scan.

## Step 1: Creating a New EVA Document and Importing Bchips.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import a Scan File dialog box will be displayed.
3. Search the Tutorials/EVA<sup>1</sup> directory and select the Bchips.RAW file.
4. Click **Open**. The scan Bchips will be displayed in the graphical view of the EVA document.




Fig. 1: Importing the Bchips.Raw file




Fig. 2: Bchips.RAW file imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.



## Step 2: Setting the Search Parameters

1. Make certain the scan is selected; if not, select it either in the Data tree or in the graphical view.
2. Click **Search / Match (scan)** in the Data Command panel  
— or —  
Right-click the scan, click **Tool** on the context menu and then **Search / Match (scan)** on the related submenu.

The **Search / Match (scan)** dialog box will be displayed.




Fig. 3: Search/match dialog box


3. Keep the default parameters:
  - ⇒ In the Candidate List tab:
    - Group duplicates
    - Criterion = **2 (Neutral)**
    - **Whole Range** selected
    - Auto



### NOTE


The Group duplicates option is not used in this example but this is a practical option which can be used without any risk.

⇒ In the Chemical Filter tab:



No chemical filter, except D (deuterium) is Discarded (red)

⇒ In the Database Filter tab:



- all Subfiles are selected
- all Quality Marks are selected
- all Status are selected except  **A** (alternative) that is cleared

4. Click the **Match** button. The candidates are listed in the Candidate List tab.

### Step 3: Matching Procedure

From the patterns listed in the dialog box, the user chooses the best candidates to identify the unknown, by comparing the stick patterns with the peaks of the current scan.

Search / Match (scan) Bchips.RAW #1

Rebuild Chemical    Chemical Filter #1 Database    Database Filter #1

Database PLU2009: 291440 - After Filters: 213504

Chemical Filter Database Filter Candidate List Selected Candidates

Index #	%	Source	ID	Quality	Status	I/Icor	Mineral	Name	Formula
<input checked="" type="checkbox"/>	117	PLU2009	PDF 00-036-0426	Star (*)	Primary		<input checked="" type="checkbox"/>	Dolomite	Ca Mg ( C O3 )2
<input type="checkbox"/>	91	PLU2009	PDF 01-084-2065	Star (*)	Primary	2.72	<input checked="" type="checkbox"/>	Dolomite	Ca Mg0.77 Fe0.23 ( C O3 )2
<input type="checkbox"/>	91	PLU2009	PDF 04-011-9828	Indexed	Primary	2.5	<input checked="" type="checkbox"/>	Dolomite	Ca1.08 Mg0.92 ( C O3 )2
<input type="checkbox"/>	91	PLU2009	PDF 04-011-9830	Indexed	Primary	2.51	<input checked="" type="checkbox"/>	Dolomite	Ca1.07 Mg0.93 ( C O3 )2
<input type="checkbox"/>	72	PLU2009	PDF 04-008-0789	Star (*)	Primary	2.53	<input checked="" type="checkbox"/>	dolomite	Ca Mg ( C O3 )2
<input type="checkbox"/>	71	PLU2009	PDF 04-011-9834	Indexed	Primary	2.77	<input checked="" type="checkbox"/>	Ankerite	Ca0.991 Mg0.676 Fe0.333 ( C O3 )2
<input type="checkbox"/>	71	PLU2009	PDF 01-074-7798	Star (*)	Primary	2.77	<input checked="" type="checkbox"/>	Ankerite	Ca ( Mg0.67 Fe0.33 ) ( C O3 )2
<input type="checkbox"/>	67	PLU2009	PDF 00-005-0622	Blank	Deleted		<input checked="" type="checkbox"/>	Dolomite	Ca Mg ( C O3 )2
<input type="checkbox"/>	57	PLU2009	PDF 00-035-0667	Star (*)	Primary		<input checked="" type="checkbox"/>	Minrecordite	Ca Zn ( C O3 )2
<input type="checkbox"/>	48	PLU2009	PDF 00-001-0942	Blank	Deleted		<input checked="" type="checkbox"/>	Dolomite	Ca O · Mg O · 2 C O2

Group Duplicates Matched 90522 / 213504 Candidates in 32.6 s.

Search / Match

Whole Range  Subrange Criterion: 2: Neutral

Auto Match

Search / Match Search by Name Filter Lists

Fig. 4: Candidate List with the first pattern selected

1. In the Candidate List tab, the first pattern in the list is selected: it is Dolomite. The corresponding ghost stick pattern (sticks and dotted lines) is displayed in the graphical view.
2. Compare the stick pattern with the current scan. The first pattern 00-036-0426 matches a part of the peaks of the scan.
3. Mark this pattern by selecting the corresponding check box (or by clicking the **Check** button below the table). The pattern is associated to a color and the image of its stick patterns remains when it is no longer selected (highlighted) in the Candidate List.
4. Go down in the list using the arrow buttons below the table. Compare the stick patterns to the scan peaks. The Calcite and Quartz patterns (01-089-1304 and 01-075-8322) appear to match the scan. Select the corresponding check boxes: the images of their stick patterns remain when they are no longer selected.
5. Leave the non-matching stick patterns unchecked.



Check button






Fig. 5: Bchips scan displayed in the graphical view with the three patterns 00-036-0426, 01-89-1304 and 01-075-8322

## Results

The three main phases are found with the initial EVA Search/Match run. A residual scan must be prepared to identify the last phase.

### Step 4: Preparing the Residual Scan

The “explained areas” (peaks which were found by the search/match) given by the identified phases must be removed.

1. In the Selected Candidates tab of the Search/Match (scan) dialog box, click **Residue**.
2. Click the first selected pattern in the Selection list: the areas to be removed are displayed as ghosts in the graphical view. Adjust the width of the zone to exclude around the pattern sticks with the slider if necessary.
3. Click the **Apply** button.
4. Proceed the same way for the other selected patterns. The removed zones are displayed using the ghost color.




Fig. 6: Scan with the removed zones in red (ghost color)

### Step 5: Performing a Subsequent Search on the Residual Scan

1. In the graphical view, zoom in on the first half of the diagram ( $5^\circ$  to about  $36^\circ 2\theta$ ). Working on relatively low angles for the identification of clay materials is highly recommended.
2. Keep the default parameters:
  - Group duplicates
  - Criterion = **2 (Neutral)**
  - **Whole Range** selected
  - Auto
  - No chemical filter, except D (deuterium) is Discarded (red)
  - all Subfiles are selected
  - all Quality Marks are selected
  - all Status are selected except  **A** (alternative) that is cleared
3. In the Candidate List tab, click the **Match** button. The candidates will be listed.



#### NOTE

The Group duplicates option is not used in this example but this is a practical option which can be used without any risk.

## Step 6: Matching Procedure



Check button



1. Working only on relatively low angles for the identification of clay materials is highly recommended. Thus, zoom in on the first half of the diagram ( $5^\circ$  to about  $36^\circ 2\theta$ )
2. In the Candidate List tab, the first pattern in the list is selected: it is Gadolinium Zirconium Selenium (pattern PDF 04-003-8951). The corresponding ghost stick pattern is displayed in the graphical view.
3. Compare the stick patterns to the current scan peaks. Identify the Clinochlore (pattern PDF 01-083-1365) which ranks 5. It ranks 2 selecting the Mineral database only.
4. Mark this pattern by selecting the corresponding check box (by clicking the **Check** button below the table). The pattern is associated to a color. The image of its stick patterns remains when it is no longer selected in the Candidate List.
5. Leave the non-matched stick patterns unchecked.




Fig. 7: Zoom on the first half of the diagram with the Clinochlore phase identified

## Step 7: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name. (To overwrite an .EVA file, use the same name).
4. Click **Save**.




Fig. 8: Bchips.EVA document

## Performing an Automatic Search/Match Operation



### NOTE

A reference database is required to perform a Search/Match operation.

This tutorial was prepared using PDF 4+ 2009 as a reference database. You may have another reference database.

If you have no database, you cannot perform this part of the tutorial concerning the Search/Match.

The following 4-step procedure describes the basic procedure for running EVA Search/Match using the automatic mode.

An Automatic Search is likely to deliver accurate results if the three following conditions are fulfilled:

- low overlap between phases
- every phase in the unknown shall have relative intensities matching the ones of its reference pattern in the database
- no phase showing a significant line broadening

On the contrary, the interactive Search/Match is quite insensitive to line overlaps, to relative intensity mismatches, whereas line broadening is not a major issue, which makes the identification of phases possible in complex scans, even with relatively strong preferred orientations.

One has to be very careful using the automatic search. It is great when the three conditions are fulfilled, but this can hardly be guessed before trying. In praxis, one will often have to keep on searching with the interactive Search/Match after a trial with the Automatic Search.

The scans used are held as tutorial files, m1.RAW, CPD-2.RAW and Mineral\_Mixture.BRML found in the Tutorial directory.

### Steps

---

1. Creating a new EVA document and importing the scan
  2. Setting of the search parameters and performing the initial search
  3. Automatic matching procedure:
    - Identification of the compounds
    - Checking with a "normal" matching procedure (for the m1 scan)
  4. Saving the EVA document containing the scan.
-



## Case #1: m1.RAW

M1 is a mixture of Calcite, Aragonite and Brucite.

### Step 1: Creating a New EVA Document and Importing m1.RAW



New button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.



Import from file button

2. Click the **Import from file** data command or the **Import from file** button. The Import a Scan File dialog box will be displayed.

3. Search the Tutorials/EVA<sup>1</sup> directory and select the m1.RAW file.

4. Click **Open**. The scan m1 will be displayed in the graphical view of the EVA document.




Fig. 9: Importing the m1.RAW file




Fig. 10: m1.RAW file imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Setting the Search Parameters

1. Make certain the scan is selected; if not, select it either in the Data tree or in the graphical view.
2. Click **Search / Match (scan)** in the Data Command panel  
— or —  
Right-click the scan, click **Tool** on the context menu and then **Search / Match (scan)** on the related submenu.

The **Search / Match (scan)** dialog box will be displayed.




Fig. 11: Search/match dialog box


3. Select the **Auto** check box.
4. Set the parameters to use for the Search/Match:
  - ⇒ In the Candidate List tab:
    - Group duplicates
    - Criterion = **2 (Neutral)**
    - **Whole Range** selected



### NOTE

The Group duplicates option is not used in this example but this is a practical option which can be used without any risk.

⇒ In the Chemical Filter tab:



No chemical filter, except D (deuterium) is Discarded (red)

⇒ In the Database Filter tab:

Search / Match (scan) m1.RAW #1

Rebuild Chemical + - [a]e Chemical Filter #1 Database + - [a]e Database Filter #1

Database PLU2009: 291440 - After Filters: 213504

Chemical Filter Database Filter Candidate List Selected Candidates

**1 Database**

- Crystallography Open Data
- My description
- PDF-4+ 2008 RDB
- PDF-4+ 2009 RDB
- PDF-4+ 2010 RDB

Filter	Value	Candidates
<input checked="" type="checkbox"/> Skip Non Ambient		264754
<input checked="" type="checkbox"/> Quality Marks		
<input checked="" type="checkbox"/> Status		
<input checked="" type="checkbox"/> Primary		212656
<input checked="" type="checkbox"/> Deleted		15723
<input type="checkbox"/> Alternate		63061
<input checked="" type="checkbox"/> Colors		
<input checked="" type="checkbox"/> Sources		
<input checked="" type="checkbox"/> Subfiles		
<input checked="" type="checkbox"/> Element # in Formula		
<input checked="" type="checkbox"/> Density		

- all Subfiles are selected
- all Quality Marks are selected
- all Status are selected except  **A** (alternative) that is cleared

5. Click the **Match** button. The candidates are listed in the Candidate List tab.

### Step 3: Automatic Matching Procedure




Fig. 12: Candidate List with the identified patterns selected



Pattern checked

The results of the Automatic Search/Match are displayed in the candidate list. The three patterns identified are automatically checked and therefore added to the data tree and to the graphical view.




Fig. 13: m1 scan displayed in the graphical view with the three patterns 01-071-3699, 00-044-1482 and 00-001-0628

### Results

The three phases found with the automatic EVA Search/Match run match well the scan and corresponds to the expected results. Nevertheless it can be interesting to check with a “non-automatic” search/match.

## “Normal” matching procedure

Set the same parameters than for the automatic search but clear the **Auto** check box before running the search.

Search / Match (scan) m1.RAW #1

Rebuild Chemical    Chemical Filter #1 Database    Database Filter #1

Database PLU2009: 291440 - After Filters: 213504

Chemical Filter Database Filter Candidate List Selected Candidates

Index #	FOM	%	Source	ID	Quality	Status	I/Icor	Mineral	Name	Formula
<input checked="" type="checkbox"/>	1	52.1	82	PLU2009	PDF 01-071-3699	Star (*)	Primary	3.23	Calcite, syn	Ca ( C O3 )
<input type="checkbox"/>	2	49.6	79	PLU2009	PDF 00-005-0586	Star (*)	Primary	2	Calcite, syn	Ca C O3
<input type="checkbox"/>	3	48.8	80	PLU2009	PDF 01-083-1762	Blank	Deleted	3.25	Calcite	Ca ( C O3 )
<input type="checkbox"/>	4	46.5	75	PLU2009	PDF 04-007-8659	Indexed	Primary	3.2	Calcite, syn	Ca ( C O3 )
<input type="checkbox"/>	5	42.6	66	PLU2009	PDF 00-047-1743	Calculated	Primary		Calcite	Ca C O3
<input type="checkbox"/>	6	32.8	37	PLU2009	PDF 00-024-0027	Calculated	Deleted		Calcite	Ca C O3
<input type="checkbox"/>	7	29.4	43	PLU2009	PDF 01-085-1108	Indexed	Deleted	3.39	Calcite, syn	Ca C O3
<input type="checkbox"/>	8	20.3	34	PLU2009	PDF 00-029-1344	Indexed	Primary		Carlinite, syn	Tl2 S
<input type="checkbox"/>	9	19.0	25	PLU2009	PDF 01-089-1304	Star (*)	Primary	3.12	Calcite, magnesium, syn	( Mg0.03 Ca0.97 ) ( C O3 )
<input type="checkbox"/>	10	19.0	31	PLU2009	PDF 00-001-0837	Blank	Deleted		Calcite	Ca C O3
<input type="checkbox"/>	11	18.4	25	PLU2009	PDF 04-012-6929	Star (*)	Primary	2.99	Calcium Magnesium Carb...	Ca0.845 Mg0.155 ( C O3 )
<input type="checkbox"/>	12	18.0	38	PLU2009	PDF 00-043-0943	Indexed	Primary		Strontium Germanium Oxide	Sr Ge O3
<input checked="" type="checkbox"/>	13	17.9	13	PLU2009	PDF 00-001-0628	Blank	Deleted		Aragonite	Ca C O3
<input type="checkbox"/>	14	15.1	13	PLU2009	PDF 00-002-0623	Indexed	Deleted		Calcite	Ca C O3
<input type="checkbox"/>	15	13.8	24	PLU2009	PDF 04-001-5486	Prototyping	Primary	9.53	Zirconium Vanadium Germ...	Zr0.91 V0.09 Ge2
<input type="checkbox"/>	16	13.7	14	PLU2009	PDF 00-024-0025	Low precision	Deleted		Aragonite	Ca C O3
<input type="checkbox"/>	17	13.0	22	PLU2009	PDF 00-044-0663	Star (*)	Primary	3.5	Arsenic Iodide Telluride	As5 Te7 I
<input type="checkbox"/>	18	12.3	43	PLU2009	PDF 01-082-1460	Star (*)	Primary	3.42	Rubidium Iron Selenate	Rb Fe ( Se O4 )2
<input type="checkbox"/>	19	12.2	14	PLU2009	PDF 01-071-4433	Indexed	Primary	4.79	Iron Vanadium Germanium	Fe V2 Ge
<input type="checkbox"/>	20	12.0	15	PLU2009	PDF 00-041-1475	Star (*)	Primary	1	Aragonite	Ca C O3

Group Duplicates Matched 69615 / 213504 Candidates in 20.7 s.

Search / Match

Whole Range Subrange Criterion: 2: Neutral

Auto Match

Search / Match Search by Name Filter Lists

Fig. 14: Candidate List with the first two patterns identified

From the patterns listed in the dialog box, the user chooses the best candidates to identify the unknown, by comparing the stick patterns with the peaks of the current scan.

1. In the Candidate List tab, the first pattern in the list is selected: it is Calcite. The corresponding ghost stick pattern (sticks and dotted lines) is displayed in the graphical view.
2. Compare the stick pattern with the current scan. The first pattern 01-071-3699 matches a part of the peaks of the scan.
3. Mark this pattern by selecting the corresponding check box (or by clicking the **Check** button below the table). The pattern is associated to a color and the image of its stick patterns remains when it is no longer selected (highlighted) in the Candidate List.



Check button



4. Go down in the list using the arrow buttons below the table. Compare the stick patterns to the scan peaks. The Aragonite and Brucite patterns (00-001-0628 and 04-011-5938) appear to match the scan. Select the corresponding check boxes: the images of their stick patterns remain when they are no longer selected.
5. Leave the non-matching stick patterns unchecked.

**Results**

The three phases found with the Automatic Search/Match are the same than those identified with the "normal" Search/Match.

## Case #2: CPD-2.RAW

CPD-2 is a mixture of Zincite, Fluorite, Al<sub>2</sub>O<sub>3</sub> and Brucite from the Rietveld quant round-robin.

### Step 1: Creating a New EVA Document and Importing CPD-2.RAW



New button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.



Import from file button

2. Click the **Import from file** data command or the **Import from file** button. The Import a Scan File dialog box will be displayed.

3. Search the Tutorials/EVA<sup>1</sup> directory and select the CPD-2.RAW file.

4. Click **Open**. The scan CPD-2 will be displayed in the graphical view of the EVA document.




Fig. 15: Importing the CPD-2.RAW file




Fig. 16: CPD-2.RAW file imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.



## Step 2: Setting the Search Parameters

1. Make certain the scan is selected; if not, select it either in the Data tree or in the graphical view.
2. Click **Search / Match (scan)** in the Data Command panel  
— or —  
Right-click the scan, click **Tool** on the context menu and then **Search / Match (scan)** on the related submenu.

The **Search / Match (scan)** dialog box will be displayed.




Fig. 17: Search/match dialog box

3. Select the **Auto** check box.
4. Set the parameters to use for the Search/Match:
  - ⇒ In the Candidate List tab:
    - Group duplicates
    - Criterion = **2 (Neutral)**
    - **Whole Range** selected



### NOTE

The Group duplicates option is not used in this example but this is a practical option which can be used without any risk.

⇒ In the Chemical Filter tab:

Search / Match (scan) CPD-2.RAW #1

Rebuild Chemical + - [a] [e] Chemical Filter #1 Database + - [a] [e] Database Filter #1

Database PLU2009: 291440 - After Filters: 213504

Chemical Filter Database Filter Candidate List Selected Candidates

Periodic table showing elements: H, D (Discarded), He, Li, Be, B, C, N, O, F, Ne, Na, Mg, Al, Si, P, S, Cl, Ar, K, Ca, Sc, Ti, V, Cr, Mn, Fe, Co, Ni, Cu, Zn, Ga, Ge, As, Se, Br, Kr, Rb, Sr, Y, Zr, Nb, Mo, Tc, Ru, Rh, Pd, Ag, Cd, In, Sn, Sb, Te, I, Xe, Cs, Ba, La, Hf, Ta, W, Re, Os, Ir, Pt, Au, Hg, Tl, Pb, Bi, Po, At, Rn, Fr, Ra, Ac, Ce, Pr, Nd, Pm, Sm, Eu, Gd, Tb, Dy, Ho, Er, Tm, Yb, Lu, Lanthanoids/Actinoids, Th, Pa, U, Np, Pu, Am, Cm, Bk, Cf, Es, Fm, Md, No, Lr.

No chemical filter, except D (deuterium) is Discarded (red)

⇒ In the Database Filter tab:

Search / Match (scan) CPD-2.RAW #1

Rebuild Chemical + - a/e Chemical Filter #1 Database + - a/e Database Filter #1

Database PLU2009: 291440 - After Filters: 21679

Chemical Filter Database Filter Candidate List Selected Candidates

1 Database

- Crystallography Open Data
- My description
- PDF-4+ 2009 RDB

Filter	Value	Candidates
<input checked="" type="checkbox"/> Skip Non Ambient		264754
<input checked="" type="checkbox"/> Quality Marks		
<input checked="" type="checkbox"/> Status		
<input checked="" type="checkbox"/> Colors		
<input checked="" type="checkbox"/> Sources		
<input type="checkbox"/> Subfiles		
<input type="checkbox"/> Organic		33466
<input type="checkbox"/> Common phase		17549
<input type="checkbox"/> Pharmaceutical and related materials		4647
<input type="checkbox"/> Inorganic		262365
<input type="checkbox"/> Forensic		10965
<input type="checkbox"/> Pigment dye		1257
<input type="checkbox"/> Educational pattern		1071
<input type="checkbox"/> Excipients and related material		1892
<input checked="" type="checkbox"/> Mineral		34212
<input type="checkbox"/> Metal & Alloy		104429
<input type="checkbox"/> Explosive		1063
<input type="checkbox"/> Battery subfile		774
<input type="checkbox"/> Ceramic (all)		9129
<input type="checkbox"/> Cement and Hydration Product		1183
<input type="checkbox"/> Ionic conductor		1216
<input type="checkbox"/> Superconducting material		6634
<input type="checkbox"/> Zeolite		3155
<input type="checkbox"/> Polymer		1138
<input type="checkbox"/> NBS pattern		2098
<input type="checkbox"/> Intercalate		390
<input type="checkbox"/> Material with bioactivity		144
<input checked="" type="checkbox"/> Element # in Formula		
<input checked="" type="checkbox"/> Density		

- only the Mineral subfile is selected
- all Quality Marks are selected
- all Status are selected except  **A** (alternative) that is cleared

5. Click the **Match** button. The candidates are listed in the Candidate List tab.

### Step 3: Automatic Matching Procedure





Fig. 18: Candidate List with the identified patterns selected

 The results of the Automatic Search/Match are displayed in the candidate list. The four patterns identified are automatically checked and therefore added to the data tree and to the graphical view.




Fig. 19: CPD-2 scan displayed in the graphical view with the four patterns 01-070-8070, 04-002-2191, 00-043-1484 and 01-071-5972

### Results

The four phases found with the automatic EVA Search/Match run match well the scan and correspond to the expected results.

### **Case #3: Mineral mixture**

Mineral mixture is a mixture of Zincite, Quartz, Calcite, Galena, Fluorite, Corundum and Rutil.

#### **Step 1: Creating a New EVA Document and Importing Mineral\_mixture.brml**



New button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.



Import from file button

2. Click the **Import from file** data command or the **Import from file** button. The Import a Scan File dialog box will be displayed.

3. Search the Tutorials/EVA<sup>1</sup> directory and select the Mineral\_mixture.brml file.

4. Click **Open**. The scan Mineral mixture will be displayed in the graphical view of the EVA document.




Fig. 20: Importing the Mineral\_mixture.brml file




Fig. 21: Mineral\_mixture.brml file imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Setting the Search Parameters

1. Make certain the scan is selected; if not, select it either in the Data tree or in the graphical view.
2. Click **Search / Match (scan)** in the Data Command panel

— or —

Right-click the scan, click **Tool** on the context menu and then **Search / Match (scan)** on the related submenu.

The **Search / Match (scan)** dialog box will be displayed.




Fig. 22: Search/match dialog box

3. Select the **Auto** check box.
4. Set the parameters to use for the Search/Match:
  - ⇒ In the Candidate List tab:
    - Group duplicates
    - Criterion = **2 (Neutral)**
    - **Whole Range** selected



### NOTE

The Group duplicates option is not used in this example but this is a practical option which can be used without any risk.

⇒ In the Chemical Filter tab:

Search / Match (scan) Mineral\_Mixture.brml #1

Rebuild Chemical + -  $\bar{a}/e$  Chemical Filter #1 Database + -  $\bar{a}/e$  Database Filter #1

Database PLU2009: 291440

Chemical Filter Database Filter Candidate List Selected Candidates

Periodic Table showing elements with checkboxes. The 'D' isotope of Hydrogen is highlighted in red, indicating it is discarded.

Lanthanoids/Actinoids

No chemical filter, except D (deuterium) is Discarded (red).

⇒ In the Database Filter tab:

The screenshot shows the 'Database Filter' tab in the DIFFRAC.EVA software. The window title is 'Search / Match (scan) Mineral\_Mixture.brml #2'. The interface includes a 'Rebuild' button, 'Chemical' and 'Database' dropdown menus, and a 'Database Filter #1' dropdown. The main area is divided into 'Chemical Filter' and 'Database Filter' tabs. The 'Database Filter' tab is active, showing a list of filters and their corresponding candidate counts. The 'Filter' column lists various categories, and the 'Candidates' column shows the number of items in each category. The 'Value' column is empty.

Filter	Value	Candidates
<input checked="" type="checkbox"/> Skip Non Ambient		264754
<input checked="" type="checkbox"/> Quality Marks		
<input checked="" type="checkbox"/> Status		
<input checked="" type="checkbox"/> Primary		212656
<input checked="" type="checkbox"/> Deleted		15723
<input checked="" type="checkbox"/> Alternate		63061
<input checked="" type="checkbox"/> Colors		
<input checked="" type="checkbox"/> Sources		
<input type="checkbox"/> Subfiles		
<input type="checkbox"/> Organic		33466
<input type="checkbox"/> Common phase		17549
<input type="checkbox"/> Pharmaceutical and related materials		4647
<input type="checkbox"/> Inorganic		262365
<input type="checkbox"/> Forensic		10965
<input type="checkbox"/> Pigment dye		1257
<input type="checkbox"/> Educational pattern		1071
<input type="checkbox"/> Excipients and related material		1892
<input checked="" type="checkbox"/> Mineral		34212
<input type="checkbox"/> Metal & Alloy		104429
<input type="checkbox"/> Explosive		1063
<input type="checkbox"/> Battery subfile		774
<input type="checkbox"/> Ceramic (all)		9129
<input type="checkbox"/> Cement and Hydration Product		1183
<input type="checkbox"/> Ionic conductor		1216
<input type="checkbox"/> Superconducting material		6634
<input type="checkbox"/> Zeolite		3155
<input type="checkbox"/> Polymer		1138
<input type="checkbox"/> NBS pattern		2098
<input type="checkbox"/> Intercalate		390
<input type="checkbox"/> Material with bioactivity		144
<input checked="" type="checkbox"/> Element # in Formula		
<input checked="" type="checkbox"/> Density		

- only the **Mineral** subfile is selected
- all Quality Marks are selected
- all Status are selected

5. Click the **Match** button. The candidates are listed in the Candidate List tab.



### Step 3: Automatic Matching Procedure




Fig. 23: Candidate List with the found patterns checked

■ Pattern checked

The results of the Automatic Search/Match are displayed in the candidate list. The seven patterns identified are automatically checked and therefore added to the data tree and to the graphical view.




Fig. 24: Mineral\_mixture scan displayed in the graphical view with the seven patterns 03-065-0466, 01-089-7102, 01-083-1762, 00-005-0592, 01-088-2301, 04-006-2654 and 04-007-9625

### Results

The seven phases found with the automatic EVA Search/Match run match well the scan and correspond to the expected phases.

## Step 4: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name. (To overwrite an .EVA file, use the same name).
4. Click **Save**.




Fig. 25: Auto\_search.EVA document

## Creating and Using Filter Lists

The following procedure describes how to create and use filter lists.

The document used is held as a tutorial file, BX100.RAW, found in the Tutorial directory.

Steps:

---

1. Creating a new EVA document and importing BX100.RAW.
  2. Creating a filter list.
  3. Using a filter list.
- 



### NOTE

A reference database is required to perform a Search/Match operation.

This tutorial was prepared using PDF 4+ 2009 as a reference database. You may have another reference database.

If you have no database, you cannot perform this part of the tutorial using the Search/Match.

---

## Step 1: Creating a New EVA Document and Importing BX100.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed. Search the Tutorials/EVA<sup>1</sup> directory and select the BX100.RAW file.
3. Click **Open**. The scan BX100 will be displayed in the graphical view of the EVA document.




Fig. 26: Import a Scan File dialog box




Fig. 27: BX100.RAW imported in the graphical view


<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Creating a Filter List

When creating a filter list, there are three ways to add pattern to this list. Each way will be illustrated here.

### Adding the document patterns:

1. Proceed as described in the tutorial chapter “Performing the Semi-Quantitative Phase-Analysis and Comparing” on page 98 until the first Search/match step: Boehmite, Iron Oxide and Kaolinite-1A patterns identified.
2. Close the Search/Match dialog box. The three patterns are listed in the Data tree and displayed in the graphical view.
3. Select the BX100.RAW scan in the Data tree.
4. Click **Search / Match (scan)** in the Data Command panel  
— or —  
Right-click the scan, click **Tool** on the context menu and then **Search / Match (scan)** on the related submenu.  
The **Search / Match (scan)** dialog box will be displayed.
5. Go to the **Filter Lists** tab.



6. Click the **Doc. Patterns** button: the patterns from the current document are added to the Filter list.
7. Enter a name in the **Name** field (limited to 13 characters, spaces are replaced by “\_”): *LIST\_DP* (for list doc patterns).
8. Enter a description in the **Description** field (optional): *Doc Patterns BX100*.

### **i** NOTE

If no description is entered, the name will be used as the list name. If a description is entered, the description will be used as the list name in DIFFRAC.EVA (Database tab of the Settings dialog box, Database filter tab of the Search/match dialog box...).

The screenshot shows two windows from the DIFFRAC.EVA software. The top window, 'Data Tree', displays a hierarchical view of the data. Under 'Views', '2Theta View' is selected, showing 'BX100 (Coupled TwoTheta/Theta)'. Below it, 'Settings' shows '1 Chemical Filter - 1 Database Filter'. The '2Theta' section is expanded to show '1 Scan', which is further expanded to show 'Pattern List #7' containing 3 patterns: 'PDF 04-010-5683 Bohmite', 'PDF 01-077-9927 Hematite, syn', and 'PDF 00-058-2005 Kaolinite-1A'. A curved arrow points from the 'Pattern List #7' section to the 'Search / Match (scan) BX100.RAW #1' window below.

The 'Search / Match (scan) BX100.RAW #1' window shows search parameters: 'Chemical Filter #1' and 'Database Filter #1'. The database is 'PLU2009: 291440 - After Filters: 213504'. The 'Candidate List' tab is active, displaying a table of candidates:

Index #	%	Source	ID	Quality	Status	I/Icor	Mineral	Name	Formula
1	160	PLU2009	PDF 01-083-2384	Indexed	Deleted	2.72	<input checked="" type="checkbox"/>	Bohmite, syn	Al O ( O H )
2	39	PLU2009	PDF 01-089-0599	Star (*)	Alternate	3.09	<input checked="" type="checkbox"/>	Hematite, syn	Fe2 O3
3	33	PLU2009	PDF 01-080-0885	Blank	Alternate	1.05	<input checked="" type="checkbox"/>	Kaolinite-1A	Al2 ( Si2 O5 ) ( O H )4


Below the table, there are navigation buttons and a 'Listed 3 Candidates' indicator. The 'Filter Lists' section at the bottom shows a 'Filter File' named 'LIST\_DP' with a description 'Doc Patterns BX100'. The 'Subset' and 'Exclude' checkboxes are both checked. The 'Save' button is visible.

Fig. 28: Adding the document patterns to the filter list

9. Define the list type. It is recommended to select both types. Select the **Subset** and **Exclude** check boxes.
10. Click the **Save** button to save the filter list.

**Adding the Search by Name results:**

1. Select the BX100.RAW scan in the Data tree.
2. Click **Search / Match (scan)** in the Data Command panel  
— or —  
Right-click the scan, click **Tool** on the context menu and then **Search / Match (scan)** on the related submenu.  
The **Search / Match (scan)** dialog box will be displayed.
3. Go to the **Search by Name** tab.
4. In the **Names** field, enter Anatase.



5. Click the **List** button to run the search: 17 candidates are found.
6. Once the search is completed, go to the **Filter Lists** tab. Click the **Search by Name** button: the patterns from the current Search by Name results are added to the Filter list.
7. Enter a name in the **Name** field (limited to 13 characters, spaces are replaced by “\_”): *LIST\_SBN* (for list Search by Name).
8. Enter a description in the **Description** field (optional): *Search by Name Anatase*.

**i** **NOTE**

If no description is entered, the name will be used as the list name. If a description is entered, the description will be used as the list name in DIFFRAC.EVA (Database tab of the Settings dialog box, Database filter tab of the Search/match dialog box...).

9. Define the list type. It is recommended to select both types. Select the **Subset** and **Exclude** check boxes.
10. Click the **Save** button to save the filter list.

The image shows two screenshots of the DIFFRAC.EVA software interface. The top screenshot shows the 'Search / Match (scan) BX100.RAW #1' window. It displays a table of search results for 'Anatase' with columns for Index #, %, Source, ID, Quality, Status, I/Icor, Mineral, Name, and Formula. The 'List' button is visible at the bottom right of the search results section.

The bottom screenshot shows the same window with the 'Filter Lists' section expanded. The 'Filter Lists' section contains a 'List Operations' panel with 'Insert' selected and 'Delete' unselected. The 'Filter File' panel shows the name 'LIST\_SBN' and the description 'Search by Name Anatase'. The 'List Operations' panel also includes buttons for 'Selection', 'Search by Name', and 'Doc. Patterns'. The 'Filter File' panel includes buttons for 'Delete', 'Read', and 'Save'.


Index #	%	Source	ID	Quality	Status	I/Icor	Mineral	Name	Formula
1	16	PLU2009	PDF 00-001-0562	Blank	Deleted		<input checked="" type="checkbox"/>	Anatase	Ti O2
2	6	PLU2009	PDF 00-002-0387	Blank	Deleted		<input checked="" type="checkbox"/>	Anatase	Ti O2
3	10	PLU2009	PDF 00-002-0406	Low precision	Deleted		<input checked="" type="checkbox"/>	Anatase	Ti O2
4	23	PLU2009	PDF 00-004-0477	Indexed	Deleted		<input checked="" type="checkbox"/>	Anatase, syn	Ti O2
5	22	PLU2009	PDF 00-021-1272	Star (*)	Primary	3.3	<input checked="" type="checkbox"/>	Anatase, syn	Ti O2
6	7	PLU2009	PDF 00-046-1391	Indexed	Primary		<input checked="" type="checkbox"/>	Anatase, niobian	( Ti , Nb , Fe ) O2
7	13	PLU2009	PDF 01-075-2545	Star (*)	Primary	5.01	<input checked="" type="checkbox"/>	Anatase	Ti O2
8	23	PLU2009	PDF 01-075-8898	Blank	Primary	5.13	<input checked="" type="checkbox"/>	Anatase, syn	( Ti0.99 Ce0.01 ) O2
9	22	PLU2009	PDF 01-075-8899	Blank	Primary	5.22	<input checked="" type="checkbox"/>	Anatase, syn	( Ti0.98 Ce0.02 ) O2
10	21	PLU2009	PDF 01-075-8900	Blank	Primary	5.39	<input checked="" type="checkbox"/>	Anatase, syn	( Ti0.96 Ce0.04 ) O2

Fig. 29: Adding the Search by Name results to the filter list




**Adding the Search/Match results:**

1. Select the BX100.RAW scan in the Data tree.
2. Click **Search / Match (scan)** in the Data Command panel  
— or —  
Right-click the scan, click **Tool** on the context menu and then **Search / Match (scan)** on the related submenu.  
The **Search / Match (scan)** dialog box will be displayed.
3. In the **Chemical Filter** tab, set all the elements except Ti and O (grey) to discarded (red).



4. In the **Database Filter** tab, select only the **Mineral** subfile.



5. In the **Candidate List** tab, keep the default parameters:

- Group duplicates
- Criterion = **2 (Neutral)**
- **Whole Range** selected
- Auto

**NOTE**

The Group duplicates option is not used in this example but this is a practical option which can be used without any risk.

6. Click the **Match** button: 32 candidates are found (with PDF-4+ 2009).
7. Once the search is completed, go to the **Filter Lists** tab. Click the **Search/Match** button: the patterns from the current Search/Match results are added to the Filter list.
8. Enter a name in the **Name** field (limited to 13 characters, spaces are replaced by “\_”): *LIST\_SM* (for list Search/Match).
9. Enter a description in the **Description** field (optional): *Search Match BX100*.
10. Define the list type. It is recommended to select both types. Select the **Subset** and **Exclude** check boxes.
11. Click the **Save** button to save the filter list.

**NOTE**

If no description is entered, the name will be used as the list name. If a description is entered, the description will be used as the list name in DIFFRAC.EVA (Database tab of the Settings dialog box, Database filter tab of the Search/match dialog box...).

Search / Match (scan) BX100.RAW #1

Rebuild Chemical     Chemical Filter #1 Database     Database Filter #1

Database PLU2009: 291440 - After Filters: 51

Chemical Filter Database Filter Candidate List Selected Candidates

Index #	%	Source	ID	Quality	Status	I/Icor	Mineral	Name	Formula	FOM	Mtc	nM	
<input type="checkbox"/>	1	24	PLU2009	PDF 01-086-1157	Indexed	Primary	3.6	<input checked="" type="checkbox"/>	Anatase, syn	Ti0.72 O2	4.9	10	1
<input type="checkbox"/>	2	23	PLU2009	PDF 00-004-0477	Indexed	Deleted		<input checked="" type="checkbox"/>	Anatase, syn	Ti O2	4.1	12	0
<input type="checkbox"/>	3	28	PLU2009	PDF 03-065-5714	Star (*)	Primary	5.04	<input checked="" type="checkbox"/>	Anatase, syn	Ti O2	4.1	7	5
<input type="checkbox"/>	4	22	PLU2009	PDF 00-021-1272	Star (*)	Primary	3.3	<input checked="" type="checkbox"/>	Anatase, syn	Ti O2	3.8	10	2
<input type="checkbox"/>	5	16	PLU2009	PDF 00-001-0562	Blank	Deleted		<input checked="" type="checkbox"/>	Anatase	Ti O2	3.7	8	1
<input type="checkbox"/>	6	13	PLU2009	PDF 01-075-2545	Star (*)	Primary	5.01	<input checked="" type="checkbox"/>	Anatase	Ti O2	2.9	10	1
<input type="checkbox"/>	7	6	PLU2009	PDF 01-089-0555	Indexed	Primary	3.2	<input checked="" type="checkbox"/>	Rutile, syn	Ti0.924 O2	2.9	11	1
<input type="checkbox"/>	8	6	PLU2009	PDF 01-089-0554	Indexed	Primary	3.26	<input checked="" type="checkbox"/>	Rutile, syn	Ti0.936 O2	2.7	10	2
<input type="checkbox"/>	9	7	PLU2009	PDF 01-089-0552	Indexed	Primary	3.2	<input checked="" type="checkbox"/>	Rutile, syn	Ti0.928 O2	2.2	10	2
<input type="checkbox"/>	10	15	PLU2009	PDF 00-001-1292	Blank	Deleted		<input checked="" type="checkbox"/>	Rutile	Ti O2	2.1	10	1
<input type="checkbox"/>	11	6	PLU2009	PDF 00-004-0551	Indexed	Deleted		<input checked="" type="checkbox"/>	Rutile	Ti O2	2.1	9	4
<input type="checkbox"/>	12	6	PLU2009	PDF 01-073-1782	Indexed	Primary	3.45	<input checked="" type="checkbox"/>	Rutile, syn	Ti O1.95	2.0	10	2
<input type="checkbox"/>	13	5	PLU2009	PDF 00-021-1276	Star (*)	Primary	3.4	<input checked="" type="checkbox"/>	Rutile, syn	Ti O2	2.0	11	3
<input type="checkbox"/>	14	3	PLU2009	PDF 01-086-1155	Indexed	Primary	4.05	<input checked="" type="checkbox"/>	Anatase, syn	Ti0.784 O2	1.9	10	3
<input type="checkbox"/>	15	6	PLU2009	PDF 01-086-0148	Star (*)	Primary	3.57	<input checked="" type="checkbox"/>	Rutile, syn	Ti0.992 O2	1.9	9	3

Group Duplicates       Matched 32 / 51 Candidates in 68.1 s.

Search / Match

Whole Range  Subrange Criterion: 2: Neutral

Auto

Search / Match Search by Name Filter Lists

---

Search / Match (scan) BX100.RAW #1

Rebuild Chemical     Chemical Filter #1 Database     Database Filter #1

Database PLU2009: 291440 - After Filters: 51

Chemical Filter Database Filter Candidate List Selected Candidates

Index #	%	Source	ID	Quality	Status	I/Icor	Mineral	Name	Formula	
<input type="checkbox"/>	1	24	PLU2009	PDF 01-086-1157	Indexed	Primary	3.6	<input checked="" type="checkbox"/>	Anatase, syn	Ti0.72 O2
<input type="checkbox"/>	2	23	PLU2009	PDF 00-004-0477	Indexed	Deleted		<input checked="" type="checkbox"/>	Anatase, syn	Ti O2
<input type="checkbox"/>	3	28	PLU2009	PDF 03-065-5714	Star (*)	Primary	5.04	<input checked="" type="checkbox"/>	Anatase, syn	Ti O2
<input type="checkbox"/>	4	22	PLU2009	PDF 00-021-1272	Star (*)	Primary	3.3	<input checked="" type="checkbox"/>	Anatase, syn	Ti O2
<input type="checkbox"/>	5	16	PLU2009	PDF 00-001-0562	Blank	Deleted		<input checked="" type="checkbox"/>	Anatase	Ti O2
<input type="checkbox"/>	6	13	PLU2009	PDF 01-075-2545	Star (*)	Primary	5.01	<input checked="" type="checkbox"/>	Anatase	Ti O2
<input type="checkbox"/>	7	6	PLU2009	PDF 01-089-0555	Indexed	Primary	3.2	<input checked="" type="checkbox"/>	Rutile, syn	Ti0.924 O2
<input type="checkbox"/>	8	6	PLU2009	PDF 01-089-0554	Indexed	Primary	3.26	<input checked="" type="checkbox"/>	Rutile, syn	Ti0.936 O2
<input type="checkbox"/>	9	7	PLU2009	PDF 01-089-0552	Indexed	Primary	3.2	<input checked="" type="checkbox"/>	Rutile, syn	Ti0.928 O2
<input type="checkbox"/>	10	15	PLU2009	PDF 00-001-1292	Blank	Deleted		<input checked="" type="checkbox"/>	Rutile	Ti O2
<input type="checkbox"/>	11	6	PLU2009	PDF 00-004-0551	Indexed	Deleted		<input checked="" type="checkbox"/>	Rutile	Ti O2
<input type="checkbox"/>	12	6	PLU2009	PDF 01-073-1782	Indexed	Primary	3.45	<input checked="" type="checkbox"/>	Rutile, syn	Ti O1.95
<input type="checkbox"/>	13	5	PLU2009	PDF 00-021-1276	Star (*)	Primary	3.4	<input checked="" type="checkbox"/>	Rutile, syn	Ti O2
<input type="checkbox"/>	14	3	PLU2009	PDF 01-086-1155	Indexed	Primary	4.05	<input checked="" type="checkbox"/>	Anatase, syn	Ti0.784 O2
<input type="checkbox"/>	15	6	PLU2009	PDF 01-086-0148	Star (*)	Primary	3.57	<input checked="" type="checkbox"/>	Rutile, syn	Ti0.992 O2

Group Duplicates       Listed 32 Candidates

Filter Lists

Insert  Delete

Name: LIST\_SM

Description: Search Match BX100

Subset  Exclude


Search / Match Search by Name Filter Lists

Fig. 30: Adding the Search by Name results to the filter list

### Step 3: Using a Filter List

You will be given here an example of Search/Match procedure using a filter list. Here we are looking for a specific phase which is Anatase.


1. Select the BX100.RAW scan in the Data tree.
2. Click **Search / Match (scan)** in the Data Command panel — or —  
Right-click the scan, click **Tool** on the context menu and then **Search / Match (scan)** on the related submenu.  
The **Search / Match (scan)** dialog box will be displayed.
3. In the Chemical Filter tab, keep the default parameters:



No chemical filter, except D (deuterium) is Discarded (red).

## 4. In the Database Filter tab:

- all Subfiles are selected
- all Quality Marks are selected
- all Status are selected except  **A** (alternative) that is cleared
- In addition to the PDF-4+ database, select the **Search by Name Anatase** subset list.



## 5. In the Candidate List tab, keep the default parameters:

- Group duplicates
- Criterion = **2 (Neutral)**
- **Whole Range** selected
- Auto

**NOTE**

The Group duplicates option is not used in this example but this is a practical option which can be used without any risk.

6. Click the **Match** button. 15 Anatase candidate patterns are listed (with PDF-4+ 2009).

Search / Match (scan) BX100.RAW #1

Rebuild Chemical   Chemical Filter #1 Database   Database Filter #1

Database PLU2009: 291440 - After Filters: 17

Chemical Filter Database Filter Candidate List Selected Candidates

Index #	%	Source	ID	Quality	Status	I/Icor	Mineral	Name	Formula	FOM	Mtc	nM	
<input type="checkbox"/>	1	24	PLU2009	PDF 01-086-1157	Indexed	Primary	3.6	<input checked="" type="checkbox"/>	Anatase, syn	Ti0.72 O2	4.9	10	1
<input type="checkbox"/>	2	23	PLU2009	PDF 01-075-8898	Blank	Primary	5.13	<input checked="" type="checkbox"/>	Anatase, syn	(Ti0.99 Ce0.01) O2	4.5	12	0
<input type="checkbox"/>	3	19	PLU2009	PDF 04-006-9241	Prototyping	Primary	4.89	<input checked="" type="checkbox"/>	Anatase, syn	Ti0.99 V0.01 O2	4.2	8	3
<input type="checkbox"/>	4	23	PLU2009	PDF 00-004-0477	Indexed	Deleted		<input checked="" type="checkbox"/>	Anatase, syn	Ti O2	4.1	12	0
<input type="checkbox"/>	5	28	PLU2009	PDF 03-065-5714	Star (*)	Primary	5.04	<input checked="" type="checkbox"/>	Anatase, syn	Ti O2	4.1	7	5
<input type="checkbox"/>	6	22	PLU2009	PDF 00-021-1272	Star (*)	Primary	3.3	<input checked="" type="checkbox"/>	Anatase, syn	Ti O2	3.8	10	2
<input type="checkbox"/>	7	16	PLU2009	PDF 00-001-0562	Blank	Deleted		<input checked="" type="checkbox"/>	Anatase	Ti O2	3.7	8	1
<input type="checkbox"/>	8	22	PLU2009	PDF 01-075-8899	Blank	Primary	5.22	<input checked="" type="checkbox"/>	Anatase, syn	(Ti0.98 Ce0.02) O2	3.2	11	1
<input type="checkbox"/>	9	21	PLU2009	PDF 01-075-8900	Blank	Primary	5.39	<input checked="" type="checkbox"/>	Anatase, syn	(Ti0.96 Ce0.04) O2	3.1	11	1
<input type="checkbox"/>	10	13	PLU2009	PDF 01-075-2545	Star (*)	Primary	5.01	<input checked="" type="checkbox"/>	Anatase	Ti O2	2.9	10	1

Group Duplicates       Matched 15 / 17 Candidates in 52.1 s.

Search / Match

Whole Range  Subrange Criterion: 2: Neutral  Auto

Search / Match Search by Name Filter Lists

Fig. 31: Search/Match results using a Subset list

## Creating a User Database

It is possible to create and maintain a separate database containing the user own patterns. This database is called a user database.

This chapter describes the basic procedure to create a user database and to use it when performing a search.

The used scans are m1.RAW and LaB6.BRML, stored in the Tutorial directory.

Steps:

---

1. Creating a new EVA document and importing m1.RAW.
  2. Automatic matching procedure:  
Identification of the compounds
  3. Performing the Tune cell operation on the Calcite pattern.
  4. Creating the user database.
  5. Adding the Calcite pattern to the database.
  6. Importing LaB6.BRML.
  7. Performing a peak search on LaB6 and creating the corresponding DIF.
  8. Adding the LaB6 DIF to the user database.
  9. Using the User database during a Search/Match operation.
  10. Saving the EVA document.
- 



### NOTE

A reference database is required to perform a Search/Match operation.

This tutorial was prepared using PDF 4+ 2009 as a reference database. The user may have another reference database.

If the user does not have a database, this part of the tutorial using the Search/Match cannot be performed.

---

## Step 1: Creating a New EVA Document and Importing m1.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated toolbar button. The EVA document will be empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import a Scan File dialog box will be displayed.
3. Search the Tutorials/EVA<sup>1</sup> directory and select the m1.RAW file.
4. Click **Open**. The scan m1 will be displayed in the graphical view of the EVA document.




Fig. 32: Importing the m1.RAW file




Fig. 33: m1.RAW file imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7.

To be able to access the folder via Explorer, select "Display hidden files and folders" in the Folders options.



## Step 2: Performing the Automatic Matching Procedure

1. Make sure that the scan has been selected. If the scan has not been chosen select it either in the Data tree or in the graphical view.
2. Click **Search / Match (scan)** in the Data Command panel  
— or —  
right-click the scan, click **Tool** on the context menu and then **Search / Match (scan)** on the related submenu.  
The **Search / Match (scan)** dialog box will be displayed.




Fig. 34: Search/match dialog box

3. Select the **Auto** check box.
4. Set the parameters for the Search/Match procedure:
  - ⇒ In the Candidate List tab:
    - Group duplicates
    - Criterion = **2 (Neutral)**
    - **Whole Range** selected



### NOTE


The Group duplicates option is not used in this example but this is a practical option which can be used without any risk.

⇒ In the Chemical Filter tab:

The screenshot shows the 'Search / Match (scan) CPD-2.RAW #1' window. At the top, there are controls for 'Chemical' and 'Database' filters, both currently set to 'None'. Below this, the text 'Database PLU2009: 291440 - After Filters: 213504' is displayed. The main area is titled 'Chemical Filter' and contains a periodic table of elements. Each element is represented by a small button with its symbol and a small red bar below it. The element 'D' (Deuterium) is highlighted with a red bar, indicating it is the only element selected in the filter. A large empty rectangular box is present in the upper left quadrant of the periodic table area. At the bottom left, there is a button labeled 'Lanthanoids/Actinoids'.

No chemical filter, except D (deuterium) is Discarded (red)

⇒ In the Database Filter tab:



- all Subfiles are selected
- all Quality Marks are selected
- all Status are selected except  A (alternative) which has been cleared

5. Click the **Match** button. The candidates are listed in the Candidate List tab.




Fig. 35: Candidate List with the identified patterns selected

■  
Pattern checked

The results of the Automatic Search/Match are displayed in the candidate list. The three patterns identified are automatically checked and added to the data tree and to the graphical view.





Fig. 36: m1 scan displayed in the graphical view with the three patterns 01-071-3699, 00-044-1482 and 00-001-0628


### Step 3: Performing the Tune Cell Operation

1. Select the Calcite pattern in the data tree.
2. Click **Tune Cell** in the Data Command panel  
— or —  
right-click the pattern, then click **Tune Cell** on the context menu.  
The Tune cell dialog box will be displayed.



3. Do not modify anything and click the **Replace** button.

The Calcite pattern PDF 00-005-0586 will be replaced by the pattern PDF 00-005-0586 (Tune Cell) in the graphical as well as in the data tree.



## Step 4: Creating the User Database

1. Click **User Database** in the Data Command panel  
— or —  
right-click the tuned pattern and click **User Database** on the context menu.

The User Database dialog box will be displayed.




Fig. 37: User Database dialog box

2. Click the **+** button next to the **User database** field. The User Database Creation dialog box will be displayed.





Fig. 38: User database creation dialog box

3. Enter a **Short Name** and a **Long Name**: for example, *DATABASE1* and *My first user database*. The Long name is the name which will be used in DIFFRAC.EVA (e. g. in the Database tab of the Settings dialog box and in the Database filter tab of the Search/match dialog box). Please choose the name carefully. It cannot be changed after the user database has been created.
4. Enter a **Prefix** for the user patterns' name: for example, *USR*.
5. Enter the **Digit** groups to define the way the user database patterns are numbered. Enter two digits: the first digit gives the number of digits in the first digits group. The second digit gives the number of digits in the second digits group. Enter 34. The name given to the patterns will have the form *XXX-XXXX*.
6. The example field will give an example of a pattern name with a prefix and a digits group chosen.




The 'User Database Creation' dialog box contains the following fields and buttons:

- Short Name: DATABASE1
- Long Name: My first user database
- Prefix: USR ex: ABC
- Digit groups: 34 ex: 25 meaning: XX-XXXXX
- Example: USR 345-6789
- Buttons: Cancel, Create

7. Finally, click the **Create** button.

The user database will be displayed in the User database dialog box and it will be then possible to add patterns.



The 'User Database PDF 01-071-3699 (Tune Cell)' dialog box is divided into two main sections:

**Left Section (Description, Bibliography, Experimental, Cell, Peaks):**

- Name: Calcite, syn
- Formula: Ca ( C O3 )  
*Open to test the formula's validity*
- Keywords:   
*2 or more alphabetical characters*
- Color:   
*Choose one of the colors, or type a new one*
- Subset: Inorganic  Mineral

**Right Section (User Database):**

- User Database: DATABASE1 (dropdown) + (add) X (delete)
- My first user database (text box)
- Pattern ID: USR 000-0001 (dropdown) ◀ ▶
- This pattern number is available.*
- Buttons: Update ID, Delete ID
- and Document  Append ID

### Step 5: Adding the Calcite pattern to the database

Once the user database has been created, it is possible to create user patterns.

1. In the User database dialog box, as shown at the previous step, the tuned PDF calcite pattern is given by default the first pattern number available: USR 000-0001 (**Pattern ID** field).
2. Click the **Append ID** button to add the pattern to the user database.  
The User Database dialog box will be closed.
3. Close the document. All the data will be removed from the screen.



## Step 6: Creating a New Document and Importing LaB6.BRML



New button



Import from file button

1. Click **New** on the **File** menu or use the toolbar button. The EVA document will be empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed. Search the Tutorials/EVA<sup>1</sup> directory and select the LaB6.BRML file.
3. Click **Open**. The scan LaB6 will be displayed in the graphical view of the EVA document.




Fig. 39: Import a Scan File dialog box




Fig. 40: LaB6 file imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7.

To be able to access the folder via Explorer, select "Display hidden files and folders" in the Folders options.

## Step 7: Adding the LaB6 DIF to the User Database

1. Import the LaB6.BRML scan.
2. Make sure the scan has been selected. If no scan has been chosen, select it either in the data tree or in the 1D view.
3. Click **Peak Search** in the Tool list of the Data Command panel — or — right-click the scan, and then click **Tool** on the context menu. Click **Peak Search** on the Tool submenu.

The **Peak Search** dialog box will be displayed.




Fig. 41: Peak Search dialog box and graphical view showing the ghost peaks

4. The default parameters are set in the Peak Search dialog box and the ghost peaks are displayed in the graphical view. The ghost peaks will be modified if the slider is moved.
5. Select the **Remove Ka2 contribution** check box and move the slider to adjust the peaks.





Fig. 42: Performing a peak search on the LaB6 scan

6. Click the **Make DIF** button: the corresponding DIF is added to the data tree and to the graphical view.



Data	Description
Document	
Views	
2Theta View	Calcite, Aragonite, Brucite (Coupled TwoThet...
Settings	1 Chemical Filter - 1 Database Filter
2Theta	1 Scan
LaB6.brm #1	Commander Sample ID (Coupled TwoTheta/Theta)
Pattern List #12	1 Pattern
DIF (LaB6.brm)	Commander Sample ID

7. Click **User Database** in the Data Command panel  
— or —  
right-click the DIF and click **User Database** on the context menu.

The User Database dialog box will be displayed with the user database and the first pattern number available for the pattern ID (USR 000-0002).

User Database DIF (LaB6.brm)

<p>Description <b>Bibliography</b> Experimental Cell Peaks</p> <p>Name: <input type="text" value="Commander Sample ID"/></p> <p>Formula: <input type="text"/></p> <p style="text-align: center;"><i>Open to test the formula's validity</i></p> <p>Keywords: <input type="text"/></p> <p style="text-align: center;"><i>2 or more alphabetical characters</i></p> <p>Color: <input type="text"/></p> <p style="text-align: center;"><i>Choose one of the colors, or type a new one</i></p> <p>Subset: <input type="text" value="Inorganic"/> <input type="checkbox"/> Mineral</p>	<p>User Database:</p> <p><input type="text" value="DATABASE 1"/> + X</p> <p>My first user database</p> <p>Pattern ID:</p> <p><input type="text" value="USR 000-0002"/> ◀ ▶</p> <p style="text-align: center;"><i>This pattern number is available.</i></p> <p>Update ID Delete ID</p> <p><input checked="" type="checkbox"/> and Document Append ID</p>
---	---

8. Modify the **Name** and **Formula** in the Description tab on the left. Enter respectively *Lanthanum hexaboride* and *LaB6*.

User Database DIF (LaB6.brm)

Description Bibliography Experimental Cell Peaks

Name: Lanthanum hexaboride

Formula: LaB6  
*Open to test the formula's validity*

Keywords:   
*2 or more alphabetical characters*

Color:   
*Choose one of the colors, or type a new one*

Subset: Inorganic  Mineral

User Database:  
DATABASE 1  
My first user database

Pattern ID:  
USR 000-0002  
*This pattern number is available.*


Update ID Delete ID  
 and Document Append ID

9. Click the **Append ID** button to add the new pattern to the selected user database. The User Database dialog box will be closed.
10. Delete the DIF from the data tree. It will also be removed from the graphical view.

## Step 8: Using the User Database during a Search/Match operation

A user database can be used the same way as another pattern database. The procedure below will demonstrate how to use the user database when performing a search/match.

1. Select the LaB6.BRML scan.
2. Click **Search / Match (scan)** in the Data Command panel  
— or —  
right-click the current scan, click **Tool** on the menu which appears and then **Search / Match (scan)** on the related submenu.
3. In the Chemical Filter tab, keep the default parameters: no chemical filter, except D (deuterium) is discarded (red).
4. In the Database Filter tab, select *My first user database* in addition to the PDF database (here the PDF-4+ 2009) in the Database list.



5. In the Candidate List tab, keep the default parameters:
  - Group duplicates
  - Criterion = **2 (Neutral)**
  - **Whole Range** selected
  - Auto
6. Click the **Match** button.




Fig. 43: Search Results including the user pattern

The LaB6 user pattern appears in the result list like a PDF pattern and is used in the same way in the data tree and the graphical view.




Fig. 44: User pattern added to the data tree and the graphical view

### Step 9: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name. To overwrite an .EVA file, use the same name.
4. Click **Save**.




Fig. 45: Saving LaB6\_userDB.EVA document

## Performing a Peak Search

The following chapter describes the basic procedure for running a peak search.

The scan used is held as a tutorial file, Csand.RAW, found in the Tutorial directory. Csand is white coral sand from a beach in Hawaii. It consists of three types of Carbonates, Aragonite, Calcite and Magnesian calcite.

Steps:

---

1. Creating a new EVA document and importing Csand.RAW.
  2. Setting of the setting the peak search parameters.
  3. Appending the peaks.
  4. Saving the EVA document containing the scan.
-



## Step 1: Creating a New EVA Document and Importing Csand.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed.
3. Search the Tutorials/EVA<sup>1</sup> directory and select the Csand.RAW file.
4. Click **Open**. The scan Csand.RAW will be displayed in the graphical view of the EVA document.




Fig. 46: Import a Scan File dialog box




Fig. 47: Csand scan imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Setting the Peak Search Parameters

1. Make certain the scan is selected; if not, select it either in the data tree or in the 1D view.
2. Click **Peak Search** in the Tool list of the Data Command panel  
— or —

Right-click the scan, and then click **Tool** on the context menu. Click **Peak Search** on the Tool submenu.

The **Peak Search** dialog box will be displayed.




Fig. 48: Peak Search dialog box and graphical view showing the ghost peaks

3. The default parameters are set in the Peak Search dialog box and the ghost peaks are displayed in the graphical view. If you move the slider the ghost peaks are modified. The Ka2 contribution can be removed by selecting the corresponding check box but it is not necessary here.
4. If the default parameters are satisfactory: click the **Default** button to return to the default parameters.

### Step 3: Appending the Peaks

1. To validate the peaks and add them to the graphical view as well as to the data tree, click the **Append to List** button.




Fig. 49: Peak Search dialog box and graphical view showing the ghost peaks

The screenshot shows the 'Data Tree' window in the software. The tree structure is as follows:

- Data
  - Document
    - Views
      - 2Theta View** CSAND (Coupled TwoTheta/Theta)
    - Settings
      - 1 Chemical Filter - 1 Database Filter
    - 2Theta** 1 Scan
      - Csand.raw #1 CSAND (Coupled TwoTheta/Theta)
        - Peak List #6** 34 Peaks
          - Peak #1 46.4
          - Peak #2 101
          - Peak #3 2186
          - Peak #4 1090
          - Peak #5 538
          - Peak #6 1334
          - Peak #7 109
          - Peak #8 238
          - Peak #9 1160
          - Peak #10 948

Fig. 50: Peak List added in the Data Tree

2. You can edit each peak in the Peak Property table. For example you can modify their position or change the layout.

#### **i** NOTE

You can also create a DIF pattern by clicking the **Make DIF** button.

### Step 4: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name.
4. Click **Save**.




Fig. 51: Saving Csand.EVA document

## Computing $K\alpha_2$ Stripping

The following procedure describes the basic procedure for computing  $K\alpha_2$  Stripping.

The scan used is held as a tutorial file, Quartz.RAW, found in the Tutorial directory.

Steps:

---

1. Creating a new EVA document and importing Quartz.RAW.
  2. Computing the  $K\alpha_2$  stripping.
  3. Saving the EVA document containing the scan.
-

## Step 1: Creating a New EVA Document and Importing Quartz.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed.
3. Search the Tutorials/EVA<sup>1</sup> directory and select the Quartz.RAW file.
4. Click **Open**. The scan Quartz.RAW will be displayed in the graphical view of the EVA document.




Fig. 52: Import a Scan File dialog box




Fig. 53: Quartz scan imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Computing the $K\alpha_2$ Stripping

1. Make certain the scan is selected; if not, select it either in the data tree or in the 1D view.
2. Right-click the scan, and then click **Tool** on the context menu. Click **Strip  $K\alpha_2$**  on the related submenu: the **Strip  $K\alpha_2$**  dialog box will be displayed.




Fig. 54: Strip  $K\alpha_2$  dialog box and in the graphical view: ghost  $K\alpha_2$ -stripped line

The default parameters are set in the Strip  $K\alpha_2$  dialog box and the ghost line shows the  $K\alpha_2$ -stripped line in the graphical view. If you move the slider the ghost line is modified.

3. Zoom in on the peaks to see the result clearly.




Fig. 55: Zoom on the peaks ( $2\theta=38^\circ$  to  $48^\circ$ )

4. If the default parameters are satisfactory: click the **Default** button to return to the default parameters if necessary.

**5.** Then, you can:

- Click **Append** to append the K $\alpha$ 2 subtracted scan to the document. The scan is added to the scan list and its properties can be edited in the property table.
- Click **Replace** to replace the original scan with the K $\alpha$ 2 subtracted scan.




Fig. 56: Original scan replaced with the K $\alpha$ 2 subtracted scan - Zoom on the peaks ( $2\theta=38^\circ$  to  $48^\circ$ )

**Step 3: Saving**

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name.
4. Click **Save**.




Fig. 57: Saving the Quartz.EVA document



## Smoothing Scans

The following procedure describes the basic procedure for smoothing scans.

The scan used is held as a tutorial file, Quartz.RAW, found in the Tutorial directory.

Steps:

---

1. Creating a new EVA document and importing Quartz.RAW.
  2. Smoothing the scan.
  3. Saving the EVA document containing the scan.
-

## Step 1: Creating a New EVA Document and Importing Quartz.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed.
3. Search the Tutorials/EVA<sup>1</sup> directory and select the Quartz.RAW file.
4. Click **Open**. The scan Quartz.RAW will be displayed in the graphical view of the EVA document.




Fig. 58: Import a Scan File dialog box




Fig. 59: Quartz scan imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Smoothing the Scan

### Case #1: Smoothing

1. Make certain the scan is selected; if not, select it either in the data tree or in the 1D view.
2. Right-click the scan, and then click **Tool** on the context menu. Click **Smooth** on the related submenu: the **Smooth** dialog box will be displayed.




Fig. 60: Smooth dialog box and in the graphical view: ghost smoothed line

The default parameters are set in the Smooth dialog box and the ghost line shows the smoothed line in the graphical view. If you move the slider the ghost line is modified.

3. Zoom on the peaks to clearly see the result.




Fig. 61: Zoom on peaks ( $2\theta=67^\circ$  to  $69^\circ$ )

4. If the default parameters are satisfactory: click the **Default** button to return to the default parameters.

## 5. Then, you can:

- Click **Append** to append the smoothed scan to the document. The scan is added to the scan list and its properties can be edited in the property table.
- Click **Replace** to replace the original scan with the smoothed scan.




Fig. 62: Original scan replaced with the smoothed scan - Zoom on peaks ( $2\theta=67^\circ$  to  $69^\circ$ )

### **Case #2: Fourier Smoothing**

1. Make certain the scan is selected; if not, select it either in the data tree or in the 1D view.
2. Right-click the scan, and then click **Tool** on the context menu. Click **Fourier Smooth** on the related submenu: the **Fourier Smooth** dialog box will be displayed.




Fig. 63: Smooth dialog box and in the graphical view: ghost smoothed line

3. Automatic adjustments are done. The ghost line shows the smoothed line in the graphical view. You can modify the cutoff graphically or by changing the limit value. You can modify the expansion by selecting another "Steps multiplied by" value.

#### 4. Zoom in on the peaks to see the result clearly.




Fig. 64: Zoom on peaks ( $2\theta=67^\circ$  to  $69^\circ$ )

Then, you can:

- Click **Append** to append the smoothed scan to the document. The scan is added to the scan list and its properties can be edited in the property table.
- Click **Replace** to replace the original scan with the smoothed scan.




Fig. 65: Original scan replaced with the smoothed scan - Zoom on peaks ( $2\theta=67^\circ$  to  $69^\circ$ )

### Step 3: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name.
4. Click **Save**.




Fig. 66: Saving the Smooth\_Quartz.EVA document

## Removing Aberrant Points

The following procedure describes the basic procedure for removing aberrant points from a scan. The scan used is held as a tutorial file, Spikes.RAW, found in the Tutorial directory.

Steps:

---

1. Creating a new EVA document and importing Spikes.RAW.
  2. Removing the aberrant points.
  3. Saving the EVA document containing the scan.
-

## Step 1: Creating a New EVA Document and Importing Spikes.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed.
3. Search the Tutorials/EVA<sup>1</sup> directory and select the Spikes.RAW file.
4. Click **Open**. The scan Spikes.RAW will be displayed in the graphical view of the EVA document.




Fig. 67: Import a Scan File dialog box




Fig. 68: Spikes scan imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.



## Step 2: Removing Aberrant Points

1. Make certain the scan is selected; if not, select it either in the data tree or in the 1D view.
2. Right-click the scan, and then click **Tool** on the context menu. Click **Aberrant** on the related submenu: the **Aberrant** dialog box will be displayed.




Fig. 69: Aberrant dialog box and graphical view displaying a ghost line with the aberrant points removed

The Aberrant dialog box gives the positions of the aberrant points and the corresponding uncorrected and corrected intensities. A ghost line shows the corrected scan.

3. Then, you can:

- Click **Append** to append the corrected scan to the document. The scan is added to the scan list and its properties can be edited in the property table.
- Click **Replace** to replace the original scan with the corrected scan.




Fig. 70: Original scan replaced with the corrected scan

### Step 3: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name.
4. Click **Save**.




Fig. 71: Saving the Spikes.EVA document

## Computing Areas

The following procedure describes the basic procedure for computing an area.

The scan used is held as a tutorial file, ET20.RAW, found in the Tutorial directory. ET20 is a mixture of Corundum and Boehmite.

Steps:

---

1. Creating a new EVA document and importing ET20.RAW.
  2. Selecting and computing an area.
  3. Saving the EVA document containing the scan.
-

## Step 1: Creating a New EVA Document and Importing ET20.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed.
3. Search the Tutorials/EVA<sup>1</sup> directory and select the ET20.RAW file.
4. Click **Open**. The scan ET20 will be displayed in the graphical view of the EVA document.




Fig. 72: Import a Scan File dialog box




Fig. 73: ET20.RAW file imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Selecting and Computing an Area

1. Make certain the scan is selected; if not, select it either in the data tree or in the 1D view.
2. Click **Create Area** in the Tool list of the Data Command panel,  
— or —

Right-click the scan, and then click **Tool** on the context menu. Click **Create Area** on the Tool submenu.

The **Create Area** dialog box will be displayed.

Fig. 74: Create Area dialog box

3. Click **Press and select an Area** to select an area with the mouse (press and hold the left mouse button with the pointer on one end of the selection, then point to the opposite end and release the button). In the example the first peak (Corundum peak) has been chosen.




Fig. 75: Selecting an area




Fig. 76: Area selected

You can modify the area by dragging any filled circle of the area bottom line.  
The results are displayed in the dialog box.

The dialog box 'Create Area ET20.RAW #1' displays the following data:

Angle (deg.)		Intensity (cps)	
Left	24.940	Left	28.3
Right	26.080	Right	28.6
Max	25.571	Max	1870.5
FWHM	0.149	Net Height	1842.0
Chord Mid.	25.569	Scherrer evaluation	
I. Breadth	0.181	Crystallite Size	605.809
Gravity C.	25.552	Use FWHM	<input checked="" type="radio"/>
Area (cps x deg.)		Use I. Breadth	<input type="radio"/>
Raw Area	365.5	K =	1
Net Area	333.1	Instr. Width =	0

Buttons: Press and select an Area, Append this Area

Fig. 77: Results of the area computation

- Once the selection is satisfactory, click **Append this Area** to add the area to an Area list in the Data tree.  
The results will be displayed in the Area Property table.

### Step 3: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name. (To overwrite an .EVA file, use the same name).
4. Click **Save**.




Fig. 78: Saving ET20.EVA document

## Adding Scans

The following procedure describes the basic procedure for adding scans.

The scans used are held as tutorial files, ET20.RAW, ET50.RAW, ET80.RAW, found in the Tutorial directory. They are mixtures of Corundum and Boehmite.

Steps:

---

1. Creating a new EVA document and importing ET20.RAW, ET50.RAW and ET80.RAW.
  2. Adding the scans.
  3. Saving the EVA document containing the scans.
-



## Step 1: Creating a New EVA Document and Importing ET20.RAW, ET50.RAW and ET80.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed. Search the Tutorials/EVA<sup>1</sup> directory and multi-select the ET20.RAW, ET50.RAW and ET80.RAW files.
3. Click **Open**. The scan ET20, ET50 and ET80 will be displayed in the graphical view of the EVA document.




Fig. 79: Import a Scan File dialog box




Fig. 80: ET20.RAW, ET50.RAW and ET80.RAW files imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Adding the Scans

1. Multi-select the scans in the data tree.
2. Click **Add** in the Tool list of the Data Command panel  
— or —

Right-click the multi-selection, and then click **Tool** on the context menu. Click **Add** on the Tool submenu.

The resulting scan will be displayed in the graphical view and added to the scan list in the data tree. The resulting scan will inherit the name of the first selected scan.




Fig. 81: Original scans and the scan resulting from the addition

### Step 3: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name. (To overwrite an .EVA file, use the same name).
4. Click **Save**.




Fig. 82: Saving ET\_addition.EVA document

## Subtracting Scans

The following procedure describes the basic procedure for subtracting scans. It can be applied on two scans only.

The scans used are held as tutorial files, ET20.RAW and ET50.RAW, found in the Tutorial directory. They are mixtures of Corundum and Boehmite.

Steps:

---

1. Creating a new EVA document and importing ET20.RAW and ET50.RAW.
  2. Subtracting the scans.
  3. Saving the EVA document containing the scans.
-

## Step 1: Creating a New EVA Document and Importing ET20.RAW and ET50.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed. Search the Tutorials/EVA<sup>1</sup> directory and multi-select the files ET20.RAW and ET50.RAW.
3. Click **Open**. The scans ET20 and ET50 will be displayed in the graphical view of the EVA document. Select the **Display Legend** check box to display the legend in the graphical view.




Fig. 83: Import a Scan File dialog box




Fig. 84: ET20.RAW and ET50.RAW files imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Subtracting the Scans

### Case #1: Subtraction result with positive values

1. Multi-select the scans in the data tree (select ET50.RAW first).
2. Click **Subtract** in the Tool list of the Data Command panel  
— or —  
Right-click (the multi-selection), and then click **Tool** on the context menu. Click **Subtract** on the Tool submenu.

The resulting scan is displayed in the graphical view and added to the scan list in the data tree. The resulting scan is given the name of the first selected scan.




Fig. 85: Original scans and the scan resulting from the subtraction

### Case #2: Subtraction result with negative values

1. Multi-select the scans in the data tree (select ET20.RAW first).
2. Click **Subtract** in the Tool list of the Data Command panel  
— or —  
Right-click (the multi-selection), and then click **Tool** on the context menu. Click **Subtract** on the Tool submenu.

The difference curve is automatically displayed in the extended view and added to the scan list in the data tree. The resulting scan is given the name of the first selected scan and is visible only in the extended view.




Fig. 86: Difference curve displayed in the extended view

3. Adjust the extended view scale by clicking on it and dragging up or down to zoom in or out.

### Step 3: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name. (To overwrite an .EVA file, use the same name).
4. Click **Save**.




Fig. 87: Saving ET\_subtraction.EVA document

## Merging Scans

The following procedure describes the basic procedure for subtracting scans. It can be applied on several scans.

The scans used are held as tutorial files, Corundum-1.RAW, Corundum-2.RAW and Corundum-3.RAW, found in the Tutorial directory.

### Steps:

---

1. Creating a new EVA document and importing Corundum-1.RAW, Corundum-2.RAW and Corundum-3.RAW.
  2. Merging the scans.
  3. Saving the EVA document containing the scans.
-



## Step 1: Creating a New EVA Document and Importing Corundum-1.RAW, Corundum-2.RAW and Corundum-3.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed. Search the Tutorials/EVA<sup>1</sup> directory and multi-select the files Corundum-1.RAW, Corundum-2.RAW and Corundum-3.RAW.
3. Click **Open**. The scans Corundum-1, Corundum-2 and Corundum-3 will be displayed in the graphical view of the EVA document.




Fig. 88: Import a Scan File dialog box




Fig. 89: Corundum-1.RAW, Corundum-2.RAW and Corundum-3 files imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Merging the Scans

1. Multi-select the scans from Corundum-1 to Corundum-3 in the data tree; the Corundum-1 is taken as the reference.
2. Click **Merge** in the Tool list of the Data Command panel  
— or —

Right-click the multi-selection, and then click **Tool** on the context menu. Click **Merge** on the Tool submenu.

The resulting scan is displayed in the graphical view and added to the scan list in the data tree. The resulting scan is given the name of the first selected scan.




Fig. 90: Original scans and the scan resulting from the merging

### Step 3: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name. (To overwrite an .EVA file, use the same name).
4. Click **Save**.




Fig. 91: Saving Merge\_corundum.EVA document

## Normalizing Scans

The following procedure describes the basic procedure for normalizing scans.

The scans used are held as tutorial files, Znox1.RAW and Znox2.RAW, found in the Tutorial directory. They are zinc minerals.

Steps:

---

1. Creating a new EVA document and importing Znox1.RAW and Znox2.RAW.
  2. Normalizing the scans.
  3. Saving the EVA document containing the scans.
-

## Step 1: Creating a New EVA Document and Importing Znox1.RAW and Znox2.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed. Search the Tutorials/EVA<sup>1</sup> directory and multi-select the Znox1.RAW and Znox2.RAW.
3. Click **Open**. The scans Znox1.RAW and Znox2.RAW will be displayed in the graphical view of the EVA document.




Fig. 92: Import a Scan File dialog box




Fig. 93: Znox1.RAW (in red) and Znox2.RAW (in black) files imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Normalizing the Scans

There are three different ways of normalizing scans: scans can be normalized to share the same maximum intensity, to share a common point or to share the maximum of a certain scan. All the scans displayed in the graphical view are normalized. Note that normalization is only available in CPS. Therefore if the display is in Counts, you have to switch to CPS to be able to normalize.

To normalize on the Znox2 scan:

1. The Y-scale unit is Counts by default. Select CPS for the Y-scale in the 2 Theta View Property table.
2. Right-click the Znox2 scan: the related context menu is displayed.
3. Point to **Normalize all visible scans** and on the sub-menu, click the **Normalize on scan Znox2.RAW** command.
4. The scans are normalized on the maximum of the Znox2 scan:




Fig. 94: Original scans and the scan resulting from normalizing

### Step 3: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name. (To overwrite an .EVA file, use the same name).
4. Click **Save**.




Fig. 95: Saving Znox.EVA document

## Computing the Crystallinity

The following procedure describes the basic procedure for computing the crystallinity of a sample.

The scan used is held as a tutorial file, m2.RAW, found in the Tutorial directory.

Steps:

---

1. Creating a new EVA document and importing m2.RAW.
  2. Computing the crystallinity.
  3. Checking the calculation.
  4. Saving the EVA document containing the scans.
-



## Step 1: Creating a New EVA Document and Importing m2.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed. Search the Tutorials/EVA<sup>1</sup> directory and select the m2.RAW file.
3. Click **Open**. The scan m2 will be displayed in the graphical view of the EVA document.




Fig. 96: Import a Scan File dialog box




Fig. 97: m2 file imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Computing the Crystallinity

1. Make certain the m2 scan is selected in the data tree.
2. Adjust the background. To do so:
  - Click **Background** in the Tool list of the Data Command panel  
— or —  
Right-click the scan in the data tree and then click **Tool** on the context menu. Click **Background** on the Tool submenu.  
The background dialog box will be displayed.
  - Adjust the background curvature as shown below.




Fig. 98: Background adjusted

3. Select the **Compute Crystallinity** check box in the Scan Property table.  
The crystallinity percentage is automatically computed and given below in the table.

Crystallinity	
Compute Crystallinity	<input checked="" type="checkbox"/> Yes
%-Crystallinity	39.1 %
%-Amorphous	60.9 %
Global Area	2878
Reduced Area	1125

### Step 3: Checking the Results

You can check the results by computing the global area “manually”.

1. Make certain the m2 scan is selected in the data tree.
2. Click **Background** in the Tool list of the Data Command panel  
— or —  
Right-click the scan in the data tree and then click **Tool** on the context menu. Click **Background** on the Tool submenu.
3. Set the curvature to the minimum of 0.01.
4. Click the **Append Background** button to create a scan from the background line: this scan will be named m2.RAW #2 by default.
5. Subtract this scan (m2.RAW #2) from the initial m2.RAW scan (m2.RAW #1). See the tutorial chapter “Subtracting Scans” on page 86 to know how to proceed. It gives you a new scan which is listed in the data tree as m2.RAW #3.

6. Select the scan resulting from the subtraction if it is not already selected. Compute the area on the whole range. To do so:

- Click **Create Area** in the Tool list of the Data Command panel  
— or —

Right-click the multi-selection, and then click **Tool** on the context menu. Click **Create Area** on the Tool submenu.

The Create Area dialog box will be displayed.

- Enter the left and right angle values (2 and 56.7) of the scan to compute the corresponding area. The value to look at is the raw area.

You can check if this value is equal to the global area value given automatically by the program which is 2878 cps x deg.

7. You can compute the Reduced area the same way but using the adjusted background. The raw area value found is equal to the reduced area value which is 1125 cps x deg. From this you can finally compute the crystallinity using the following formulas:

$$\%Amorphous = \frac{\text{Global area} - \text{Reduced area}}{\text{Global area}} \times 100$$

$$\%Crystallinity = 100 - \%Amorphous$$

### Step 4: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name. (To overwrite an .EVA file, use the same name).
4. Click **Save**.




Fig. 99: Saving m2.EVA document

## Simulating a Slit Mode

The following chapter describes the basic procedure for simulating a slit mode.

The scan used is held as a tutorial file, LaB6.BRML, found in the Tutorial directory.

LaB6 was measured using variable slits.

Steps:

---

1. Creating a new EVA document and importing LaB6.BRML.
  2. Performing the Search/Match operation.
  3. Simulating a slit mode.
  4. Saving the EVA document containing the scans.
-

## Step 1: Creating a New EVA Document and Importing LaB6.BRML



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed. Search the Tutorials/EVA<sup>1</sup> directory and select the LaB6.BRML file.
3. Click **Open**. The scan LaB6 will be displayed in the graphical view of the EVA document.




Fig. 100: Import a Scan File dialog box




Fig. 101: LaB6 file imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Performing the Search/Match Operation

1. Make certain the LaB6 scan is selected in the data tree.
2. Keep the default search parameters and run a search. Here are the results obtained:

Search / Match (scan) LaB6.brml #1

Rebuild Chemical    Chemical Filter #1 Database    Database Filter #1

Database PLU2009: 291440 - After Filters: 213504

Chemical Filter Database Filter Candidate List Selected Candidates

Index #	%	Source	ID	Quality	Status	I/Icor	Mineral	Name	Formula
<input checked="" type="checkbox"/>	111	PLU2009	PDF 00-059-0332	Star (*)	Primary		<input type="checkbox"/>	Lanthanum Boride	La B6
<input type="checkbox"/>	122	PLU2009	PDF 00-058-0290	Star (*)	Primary		<input type="checkbox"/>	Boron Lanthanum Samarium	B6 La0.9 Sm0.1
<input type="checkbox"/>	139	PLU2009	PDF 01-075-1402	Star (*)	Primary	9.55	<input type="checkbox"/>	Lanthanum Boron	La B5.784
<input type="checkbox"/>	141	PLU2009	PDF 01-075-1415	Star (*)	Primary	9.51	<input type="checkbox"/>	Lanthanum Boron	La B5.892
<input type="checkbox"/>	135	PLU2009	PDF 04-001-2105	Star (*)	Primary	9.53	<input type="checkbox"/>	Lanthanum Boride	La B6
<input type="checkbox"/>	135	PLU2009	PDF 01-075-1414	Star (*)	Primary	9.53	<input type="checkbox"/>	Lanthanum Boron	La B5.832
<input type="checkbox"/>	146	PLU2009	PDF 04-005-6122	Prototyping	Primary	10.2	<input type="checkbox"/>	Europium Ytterbium Boride	Eu0.3 Yb0.7 B6
<input type="checkbox"/>	116	PLU2009	PDF 01-070-8265	Blank	Primary	9.26	<input type="checkbox"/>	Lanthanum Boron	La0.983 B6
<input type="checkbox"/>	140	PLU2009	PDF 04-003-5116	Prototyping	Primary	9.4	<input type="checkbox"/>	Lanthanum Boride	La B5.90
<input type="checkbox"/>	135	PLU2009	PDF 01-078-2381	Blank	Deleted	9.53	<input type="checkbox"/>	Lanthanum Boron	La B5.83
<input type="checkbox"/>	131	PLU2009	PDF 04-004-5778	Prototyping	Primary	9.03	<input type="checkbox"/>	Lanthanum Boride	La0.97 B6
<input type="checkbox"/>	125	PLU2009	PDF 00-058-0291	Star (*)	Primary		<input type="checkbox"/>	Boron Lanthanum Samarium	B6 La0.8 Sm0.2
<input type="checkbox"/>	117	PLU2009	PDF 00-058-0336	Star (*)	Primary		<input type="checkbox"/>	Boron Lanthanum Samarium	B6 La0.7 Sm0.3

Group Duplicates

Matched 102402 / 213504 Candidates in 51.5 s.

Search / Match

Whole Range  Subrange

Criterion: 2: Neutral

Auto

Match

Search / Match Search by Name Filter Lists

Lanthanum Boride (PDF 00-059-0332) can easily be identified.

3. Mark this pattern by selecting the corresponding check box. The pattern is displayed in the graphical view and added to the data tree.
4. Adjust the Y-scale of the pattern.




Fig. 102: Zoom on the scan and the pattern ( $2\theta=20^\circ$  to  $60^\circ$ )

### Step 3: Simulating a Slit Mode

1. Make certain the LaB6 scan is selected in the data tree.
2. In the Scan Property Table, change the **Simul. Slit Mode** property from **Variable** to **Fixed**.
3. The intensities of the scan will be recalculated. The other scan related data, a pattern in this example, will be projected as well.




Fig. 103: Zoom on the scan and the pattern ( $2\theta=20^\circ$  to  $60^\circ$ )



#### Step 4: Saving

5. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
6. Select the appropriate disc drive and directory.
7. Type in the file name. (To overwrite an .EVA file, use the same name).
8. Click **Save**.




Fig. 104: Saving LaB6.EVA document

## Performing the Semi-Quantitative Phase-Analysis and Comparing with a Chemical Analysis

The following procedure describes the semi-quantitative phase-analysis on a scan and the comparison of the results with the results of a chemical analysis.

The document used is held as a tutorial file, BX100.RAW, found in the Tutorial directory.

Steps:

---

1. Creating a new EVA document and importing BX100.RAW.
  2. Performing the Search/Match operation.
  3. Exploiting the semi-quantitative phase-analysis.
  4. Comparing the Results with the Results of a Chemical Analysis
  5. Saving.
- 



### NOTE

A reference database is required to perform a Search/Match operation.

This tutorial was prepared using PDF 4+ 2009 as a reference database. You may have another reference database.

If you have no database, you cannot perform this part of the tutorial using the Search/Match.

---

## Step 1: Creating a New EVA Document and Importing BX100.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed. Search the Tutorials/EVA<sup>1</sup> directory and select the BX100.RAW file.
3. Click **Open**. The scan BX100 will be displayed in the graphical view of the EVA document.




Fig. 105: Import a Scan File dialog box




Fig. 106: BX100.RAW imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Performing the Search/Match Operation

Before the semi-quantitative phase analysis, all phases must be identified.

Here, we will describe briefly the search/match operation. For more details about how to proceed see the tutorial chapter “Performing a Search/Match Operation” on page 1 or the EVA Manual.

1. Keep the default search parameters and run a first search. Here are the results obtained:

Index #	FOM	Mtc	nM	%	Source	ID	Quality	Status	I/Icor	Mineral	Name	Formula
1	0.54	12	0	132	PLU2009	PDF 04-010-5683	Indexed	Primary	2.59	Yes	Bohmite	AlO(OH)
2	0.56	12	0	160	PLU2009	PDF 01-083-2384	Indexed	Deleted	2.72	Yes	Bohmite, syn	AlO(OH)
3	0.59	12	0	150	PLU2009	PDF 01-072-0359	Indexed	Primary	2.98	Yes	Bohmite	AlOOH
4	0.63	14	0	85	PLU2009	PDF 00-021-1307	Indexed	Primary		Yes	Bohmite, syn	AlO(OH)
5	0.72	8	0	39	PLU2009	PDF 00-020-0038	Low precision	Primary		No	Aluminum Oxide Bromide Hydroxide	Al <sub>2</sub> O <sub>3</sub> Br <sub>3</sub> (OH)
6	0.75	16	0	38	PLU2009	PDF 04-002-2984	Prototyping	Primary	2.62	No	Vanadium Iron Oxide	VFeO <sub>3</sub>
7	0.76	14	1	40	PLU2009	PDF 04-009-5898	Blank	Primary	2.58	No	Titanium Iron Oxide	Ti <sub>0.22</sub> Fe <sub>1.78</sub> O <sub>3</sub>
8	0.77	15	0	40	PLU2009	PDF 04-006-5322	Prototyping	Primary	2.6	No	Iron Aluminum Oxide	Fe <sub>1.78</sub> Al <sub>0.22</sub> O <sub>3</sub>
9	0.79	14	1	40	PLU2009	PDF 01-088-0434	Blank	Primary	2.43	No	Iron Tin Oxide	Fe <sub>1.727</sub> Sn <sub>0.205</sub> O <sub>3</sub>
10	0.79	15	0	41	PLU2009	PDF 04-006-6579	Indexed	Primary	3.04	No	Iron Oxide	Fe <sub>2</sub> O <sub>3</sub>

Bohmite (PDF 04-010-5683) and Iron Oxide (PDF 04-006-6579) can be easily identified. The search results show several types of iron oxide, we select pure iron oxide.




Fig. 107: Bohmite pattern in red and Iron Oxide pattern in blue

- Remove the already explained regions by preparing a residual scan. Another run on the residual scan makes it possible to identify Kaolinite-1A. We select the pattern PDF 04-010-4800 which has an  $I/I_{cor}$  value and belongs to a structure database. Structure databases have calculated  $I/I_{cor}$  values which are generally more reliable than the measured  $I/I_{cor}$  values.

Index #	%	Source	ID	Q	S	I/I <sub>cor</sub>	Mineral	Name
<input type="checkbox"/>	1	37	PLU2009 PDF 01-070-3507	I	P	8.08	Yes	Fourmarierite
<input type="checkbox"/>	2	9	PLU2009 PDF 00-058-2006	B	P		Yes	Kaolinite-1Ad
<input type="checkbox"/>	3	25	PLU2009 PDF 00-022-0266	O	P		No	Erbium Acetate
<input type="checkbox"/>	4	25	PLU2009 PDF 00-011-0243	I	P		No	Ammonium Hydrogen Phosphate Imide
<input type="checkbox"/>	5	9	PLU2009 PDF 00-058-2028	I	P		Yes	Kaolinite-1A
<input type="checkbox"/>	6	31	PLU2009 PDF 01-070-3517	B	P	8.79	Yes	Fourmarierite
<input type="checkbox"/>	7	8	PLU2009 PDF 00-058-2002	I	P		Yes	Dickite-2M1
<input type="checkbox"/>	8	8	PLU2009 PDF 00-058-2004	I	P		Yes	Kaolinite-1A
<input type="checkbox"/>	9	4	PLU2009 PDF 00-059-1232	I	P		No	((3,4,5-Trihydroxy-6-((3,4,5-trihydroxytetrahydro-2H-2-pyra
<input type="checkbox"/>	10	6	PLU2009 PDF 00-025-0546	B	P		Yes	Bementite
<input type="checkbox"/>	11	5	PLU2009 PDF 00-005-0145	O	P		No	1-Methyl-5-methylaminotetrazole
<input checked="" type="checkbox"/>	12	39	PLU2009 PDF 04-010-4800	I	P	0.98	Yes	Kaolinite-1A




Fig. 108: Boehmite pattern in red, Iron Oxide pattern in blue and Kaolinite-1A pattern in green.

3. Removing the regions explained by Kaolinite-1A would not be useful as the Kaolinite-1A peaks are very numerous. Using the results of a chemical analysis of the sample can help to complete the identification:

- Select the scan BX100.RAW either in the data tree or in the graphical view.
- Click **Import XRF Results** in the Data Command panel  
— or —  
Right-click the scan, then click **Import XRF Results** on the context menu.
- The Open an XRF results filename dialog box will be displayed. Select the BX100.txt file in the Open an XRF results filename dialog box and click **Open**.
- The results of the chemical analysis will be displayed in the element list in the data tree next to the SQD results.

Data	Description
BX100.RAW #1	BX100 (Coupled Theta/TwoTheta scan)
Element List #72	6 Elements
H [1] - Hydrogen	SQD=1.418 % - XRF=n.a.
O [8] - Oxygen	SQD=51.266 % - XRF=42.312 %
Al [13] - Aluminium	SQD=29.976 % - XRF=32.500 %
Si [14] - Silicon	SQD=8.313 % - XRF=3.920 %
Ti [22] - Titanium	SQD= - XRF=1.610 %
Fe [26] - Iron	SQD=9.028 % - XRF=18.300 %
Pattern List #71	3 Patterns
PDF 04-010-5683	Bohmite
PDF 04-006-6579	Iron Oxide
PDF 04-010-4800	Kaolinite-1A

The chemical analysis indicates that the sample contains titanium.

4. Create a chemical filter from Titanium. To do so:

- Select the Titanium in the data tree.
- Click **Green Filter data** in the Data Command panel.  
— or —  
Right-click the scan, then click **Create|Green Filter data** on the context menu.
- A new chemical filter is created: BX100.RAW #1 Green Filter.

5. Another run using this new chemical filter makes it then possible to identify Anatase (PDF 01-075-8898).

Index ...	%	Source	ID	Q	S	I/Icor	Mi...	Name	Formula	
<input type="checkbox"/>	1	12	PLU2009	PDF 00-031-0681	O	P	No	Lead Chromium Titanium Oxide	Pb Cr <sub>2.3</sub> Ti <sub>2.3</sub> O <sub>9</sub>	
<input checked="" type="checkbox"/>	2	23	PLU2009	PDF 01-075-8898	B	P	5.13	Yes	Anatase, syn	(Ti <sub>0.99</sub> Ce <sub>0.01</sub> ) O <sub>2</sub>
<input type="checkbox"/>	3	19	PLU2009	PDF 00-052-0325	O	P	No	Ammonium Titanium Phosphate	N H <sub>4</sub> Ti P <sub>4</sub> O <sub>13</sub>	
<input type="checkbox"/>	4	13	PLU2009	PDF 01-075-2545	S	P	5.01	Yes	Anatase	Ti O <sub>2</sub>
<input type="checkbox"/>	5	24	PLU2009	PDF 01-086-1157	I	P	3.6	Yes	Anatase, syn	Ti <sub>0.72</sub> O <sub>2</sub>
<input type="checkbox"/>	6	20	PLU2009	PDF 00-004-0477	I	D	Yes	Anatase, syn	Ti O <sub>2</sub>	
<input type="checkbox"/>	7	58	PLU2009	PDF 00-054-0889	S	P	No	Potassium Chromium Titanium Oxide Hydrate	K <sub>0.5</sub> Cr <sub>0.5</sub> Ti <sub>0.5</sub> O <sub>13</sub>	
<input type="checkbox"/>	8	18	PLU2009	PDF 01-076-2069	B	P	1.33	Yes	Vinogradovite	Na <sub>8</sub> (Ti <sub>7.44</sub> Nb <sub>0.5</sub> ) O <sub>24</sub>
<input type="checkbox"/>	9	4	PLU2009	PDF 00-057-1107	I	P	No	Ethylenediamine titanium phosphate fluoride	C <sub>2</sub> H <sub>10</sub> N <sub>2</sub> · Ti <sub>2</sub> (P O <sub>4</sub> ) <sub>2</sub>	
<input type="checkbox"/>	10	59	PLU2009	PDF 00-040-0162	O	P	No	Copper Titanium Hydrogen Phosphate Hydrate	Ti H Cu <sub>0.5</sub> (P O <sub>4</sub> ) <sub>2</sub>	
<input type="checkbox"/>	11	11	PLU2009	PDF 01-089-6535	B	P	1.29	No	Sodium Titanium Silicate	Na <sub>2</sub> (Ti Si <sub>5</sub> O <sub>13</sub> )
<input type="checkbox"/>	12	38	PLU2009	PDF 01-088-1385	B	P	5.79	No	Titanium Niobium Oxide Chloride	Ti <sub>2</sub> ((Nb <sub>6</sub> Cl <sub>6</sub> ) O <sub>4</sub> ) <sub>2</sub>




Fig. 109: Boehmite pattern in red, Iron Oxide pattern in blue, Kaolinite-1A pattern in green and Anatase pattern in pink.

### Step 3: Performing the Semi-Quantitative Phase Analysis

Once all the phases have been identified, the semi-quantitative analysis can be performed

1. The first phase identified in the list is Boehmite. Zoom in on the strongest peak of the Boehmite, and then change the Y-Scale of the pattern to adjust the stick height to the peak. To adjust the Y-Scale, point to the stick of interest and press the control key to change the pointer into a hand. Move it up or down to adjust the Y-Scale.




Fig. 110: Zoom on the peaks ( $2\theta=38^\circ$  to  $48^\circ$ )

2. Proceed the same way for the other phases.
3. Select the BX100 scan either in the data tree or in the 1D view.

4. The phases' concentrations are listed in the **S-Q** column of the pattern column view.

Scan	SS-VVV-PPPP	Compound Name	Formula	Y-Scale	I/Ic DB	I/Ic User	S-Q
BX100.RAW #1	PDF 04-010-5683	Bohmite	Al O ( OH )	99.98 %	2.590		64.29 %
BX100.RAW #1	PDF 04-006-6579	Iron Oxide	Fe2 O3	31.73 %	3.040		17.38 %
BX100.RAW #1	PDF 04-010-4800	Kaolinite-1A	Al2 Si2 O5 ( OH )4	9.10 %	0.980		15.47 %
BX100.RAW #1	PDF 01-075-8898	Anatase, syn	( Ti0.99 Ce0.01 ) O2	8.79 %	5.130		2.85 %

Fig. 111: Semi-quantitative analysis results in the Pattern Column view

To display the results as a pie or bar chart:

- Select the pattern list in the data tree.
- Click **Pattern Chart View** in the Create list of the Data Command panel  
— or —  
Right-click and then click **Create**. Click **Pattern Chart View** on the submenu.
- Select either **Pie Chart** or **Bar Chart** in the Pattern Chart View table.




Fig. 112: Semi-quantitative analysis results as a pie chart



### Step 4: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name. (To overwrite an .EVA file, use the same name).
4. Click **Save**.




Fig. 113: Saving BX100.EVA document

## Using the d Multiplied By Tool

The following procedure describes how to use the d Multiplied By tool.

The document used is held as a tutorial file, CSand.RAW, found in the Tutorial directory.

Csand is white coral sand from a beach in Hawaii.

The coral sand consists of three types of Carbonates, Aragonite, Calcite and Magnesian calcite. The calcite peaks are shifted to slightly larger angles, indicating minor Mg for Ca solid solution replacement. The Magnesian calcite peaks show large shifts.

Steps:

---

1. Creating a new EVA document and importing Csand.RAW.
  2. Performing the Search/Match operation.
  3. Saving.
- 



### NOTE

A reference database is required to perform a Search/Match operation.

This tutorial was prepared using PDF 4+ 2009 as a reference database. You may have another reference database.

If you have no database, you cannot perform this part of the tutorial using the Search/Match.

## Step 1: Creating a New EVA Document and Importing Csand.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed. Search the Tutorials/EVA<sup>1</sup> directory and select the Csand.RAW file.
3. Click **Open**. The scan Csand will be displayed in the graphical view of the EVA document.




Fig. 114: Import a Scan File dialog box




Fig. 115: Csand.RAW imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Performing the Search/Match Operation

1. Keep the default search parameters and run a first search. Here are the results obtained:

Search / Match (scan) Csand.RAW #1

Rebuild Chemical  Chemical Filter #1 Database  Database Filter #1

Database PLU2009: 291440 - After Filters: 213504

Chemical Filter Database Filter Candidate List Selected Candidates

Index #	FOM	Mtc	nM	%	Source	ID	Quality	Status	I/Icor	Mineral	Name	Formula
<input checked="" type="checkbox"/>	1	0.26	22	0	118	PLU2009 PDF 00-024-0025	Low precision	Deleted		Yes	Aragonite	Ca C O3
<input type="checkbox"/>	2	0.26	22	0	122	PLU2009 PDF 04-008-5421	Star (*)	Primary	1.14	Yes	aragonite	Ca ( C O3 )
<input type="checkbox"/>	3	0.27	20	0	121	PLU2009 PDF 01-071-2392	Star (*)	Primary	1.13	Yes	Aragonite	Ca ( C O3 )
<input type="checkbox"/>	4	0.27	25	0	125	PLU2009 PDF 00-041-1475	Star (*)	Primary	1	Yes	Aragonite	Ca C O3
<input type="checkbox"/>	5	0.29	21	0	115	PLU2009 PDF 01-070-9854	Star (*)	Primary	0.33	Yes	Aragonite	Ca ( C O3 )
<input type="checkbox"/>	6	0.38	18	0	61	PLU2009 PDF 00-001-0628	Blank	Deleted		Yes	Aragonite	Ca C O3
<input type="checkbox"/>	7	1.04	15	2	88	PLU2009 PDF 01-074-4669	Blank	Primary	10.72	No	Barium Indium Titanium Oxide	Ba2 ( ( In1.85 Ti0.15 ) O5.075 )
<input type="checkbox"/>	8	1.07	18	1	100	PLU2009 PDF 01-074-4668	Blank	Primary	10.81	No	Barium Indium Titanium Oxide	Ba2 ( ( In1.9 Ti0.1 ) O5.05 )
<input type="checkbox"/>	9	1.12	5	0	109	PLU2009 PDF 04-005-3806	Indexed	Primary	11.16	No	Titanium Niobium Carbide	Ti0.18 Nb0.82 C0.5
<input type="checkbox"/>	10	1.12	37	6	55	PLU2009 PDF 04-002-1441	Prototyping	Primary	1.22	No	Yttrium Molybdenum Iron	Y3.36 Mo0.888 Fe27.752

Group Duplicates

Matched 146505 / 213504 Candidates in 30.3 s.

Search / Match

Whole Range  Subrange

Criterion: 2: Neutral

Auto

Match

Search / Match Search by Name Filter Lists

Aragonite (PDF 00-024-0025) can be easily identified.




Fig. 116: Aragonite pattern in red

2. A run performed on the Mineral database makes it possible to find the Magnesian Calcite (PDF 00-043-0697).

Index #	%	Source	ID	Q	S	I/Icor	Mineral	Name	Formula	
<input checked="" type="checkbox"/>	1	118	PLU2009	PDF 00-024-0025	O	D		Yes	Aragonite	Ca C O3
<input type="checkbox"/>	2	122	PLU2009	PDF 04-008-5421	S	P	1.14	Yes	aragonite	Ca ( C O3 )
<input type="checkbox"/>	3	121	PLU2009	PDF 01-071-2392	S	P	1.13	Yes	Aragonite	Ca ( C O3 )
<input type="checkbox"/>	4	125	PLU2009	PDF 00-041-1475	S	P	1	Yes	Aragonite	Ca C O3
<input type="checkbox"/>	5	115	PLU2009	PDF 01-070-9854	S	P	0.33	Yes	Aragonite	Ca ( C O3 )
<input type="checkbox"/>	6	61	PLU2009	PDF 00-001-0628	B	D		Yes	Aragonite	Ca C O3
<input type="checkbox"/>	7	28	PLU2009	PDF 00-003-0405	B	D		Yes	Aragonite	Ca C O3
<input type="checkbox"/>	8	38	PLU2009	PDF 00-044-1090	I	P		Yes	Antimony Palladium	Pd31 Sb12
<input type="checkbox"/>	9	54	PLU2009	PDF 00-003-1067	O	D		Yes	Aragonite	Ca C O3
<input type="checkbox"/>	10	83	PLU2009	PDF 04-007-0914	B	P	7.68	Yes	Stibiopalladinite, syn	Pd5 Sb2
<input type="checkbox"/>	11	13	PLU2009	PDF 00-046-1430	B	P		Yes	Watanabeite	Cu4 ( As , Sb )2 S5
<input checked="" type="checkbox"/>	12	48	PLU2009	PDF 00-043-0697	S	P		Yes	Calcite, magnesian	( Ca , Mg ) C O3




Fig. 117: Aragonite pattern in red and Magnesian calcite in blue

### Step 3: Using the d Multiplied By Tool

1. Insert the pure Calcite pattern (PDF 00-005-0586) by using the Search by name tool.
2. Zoom in around the location of the 100% line of pure Calcite.
3. Select the Calcite pattern in the data tree.
4. Click **d x by** in the Data Command panel  
— or —  
Right-click the pattern, then click **d x by** on the context menu.  
The **d x by** dialog box will be displayed.
5. Set the d Multiplied by value to 0.990 using the slider. It shows the Magnesian calcite pattern (PDF 00-043-0697) corresponds to the pure Calcite d-spacing multiplied by about 0.990.





Fig. 118: Adjusting the d Multiplied by value to 0.990 using the slider

But the file shows a mixture of two solid solutions exactly on both sides. Thus, the perfect explanation of this unknown scan requires two modified Calcite patterns. We can, for instance, make use of the pure Calcite (pattern 00-005-586) and modify it by means of the **d x By** tool. The adjustment of the slider leads to a  $d$  times value of 0.984 for matching the modified Calcite pattern having the stronger peaks, and of 0.996 for the second pattern as shown in the figure below.



1 Calcite pattern 00-005-0586

2 Calcite pattern 00-005-0586  
multiplied by 0.996.

3 Magnesian Calcite pattern 00-043-0687

4 Calcite pattern 00-005-0586  
multiplied by 0.984.

### Step 4: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name. (To overwrite an .EVA file, use the same name).
4. Click **Save**.




Fig. 119: Saving Csand\_dxBy.EVA document

## Using the Tune Cell Tool

The following procedure describes how use the Tune Cell tool.

The document used is held as a tutorial file, Francolite.RAW, found in the Tutorial directory.

The measured sample is a Carbonate-fluoroapatite (its mineral name is Francolite) and the difference from the Fluorapatite  $\text{Ca}_{10}(\text{PO}_4)_6\text{F}_2$  is that P is partly replaced by C.

The replacement of  $\text{P}^{5+}$  (ionic radius 0.35) by the smaller  $\text{C}^{4+}$  (0.16) is expressed by deformation of the  $\text{PO}_4$ -tetraeder, with shorter P-O-distances and the reduction of  $a_0$  cell parameter<sup>1</sup>.

Steps:

---

1. Creating a new EVA document and importing Francolite.RAW.
  2. Performing the Search/Match operation.
  3. Performing the Tune cell operation.
  4. Saving.
- 



### NOTE

A reference database is required to perform a Search/Match operation.

This tutorial was prepared using PDF 4+ 2009 as a reference database. You may have another reference database.

If you have no database, you cannot perform this part of the tutorial using the Search/Match.

---

<sup>1</sup> X-Ray Powder Diffraction Study of Francolite by the Rietveld Method, B. Perdikatsis, Materials Science Forum Vols. 79-82 (1991), pp. 809-814.



## Step 1: Creating a New EVA Document and Importing Francolite.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed. Search the Tutorials/EVA<sup>1</sup> directory and select the Francolite.RAW file.
3. Click **Open**. The scan Francolite will be displayed in the graphical view of the EVA document.




Fig. 120: Import a Scan File dialog box




Fig. 121: Francolite.RAW imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Performing the Search/Match Operation

1. Keep the default search parameters and run a search. Here are the results obtained:

Database PLU2009: 291440 - After Filters: 213504

Chemical Filter Database Filter Candidate List Selected Candidates

Index #	%	Source	ID	Q	S	I/Icor	Mineral	Name	Formula	
<input type="checkbox"/>	1	97	PLU2009	PDF 00-031-0267	B	P		Yes	Carbonatefluorapatite	Ca <sub>10</sub> (P O <sub>4</sub> ) <sub>5</sub> C O <sub>3</sub> F <sub>1.5</sub> (O H) <sub>0.5</sub>
<input type="checkbox"/>	2	97	PLU2009	PDF 00-021-0141	B	D		Yes	Carbonate-fluorapatite	Ca <sub>10</sub> (P O <sub>4</sub> ) <sub>5</sub> C O <sub>3</sub> F <sub>1.5</sub> (O H) <sub>0.5</sub>
<input type="checkbox"/>	3	18	PLU2009	PDF 00-006-0454	B	D		Yes	Fluorellestadite, phosphatian	Ca <sub>5</sub> ((P, S, Si) O <sub>4</sub> ) <sub>3</sub> (O H, C O <sub>3</sub> )
<input checked="" type="checkbox"/>	4	37	PLU2009	PDF 00-015-0876	S	P	1.5	Yes	Fluorapatite, syn	Ca <sub>5</sub> (P O <sub>4</sub> ) <sub>3</sub> F
<input type="checkbox"/>	5	29	PLU2009	PDF 00-056-1014	O	P		No	Calcium Europium Fluoride Phosphate	Ca <sub>5</sub> (P O <sub>4</sub> ) <sub>3</sub> F - Eu
<input type="checkbox"/>	6	54	PLU2009	PDF 00-047-0262	O	P		No	Sodium Calcium Hydrogen Carbonate ...	Ca <sub>8</sub> H <sub>2</sub> (P O <sub>4</sub> ) <sub>6</sub> · H <sub>2</sub> O - Na H C O <sub>3</sub> - H <sub>2</sub> O
<input type="checkbox"/>	7	91	PLU2009	PDF 00-053-0231	S	P	4.6	No	Strontium Gallium Manganese Oxide	Sr <sub>2</sub> Ga Mn O <sub>4</sub> 97
<input type="checkbox"/>	8	36	PLU2009	PDF 00-032-0163	S	P		No	Calcium Gadolinium Oxide Phosphate	Ca <sub>8</sub> Gd <sub>2</sub> (P O <sub>4</sub> ) <sub>6</sub> O <sub>2</sub>
<input type="checkbox"/>	9	133	PLU2009	PDF 00-020-0506	O	P		No	Ammonium Iron Phosphate Hydrate	(N H <sub>4</sub> ) <sub>2</sub> Fe P <sub>2</sub> O <sub>7</sub> · 2H <sub>2</sub> O

Group Duplicates  Matched 114935 / 213504 Candidates in 19.9 s.

Carbonatefluorapatite is easily identified but we select Fluorapatite (PDF 00-015-0876) for the need of the example.





Fig. 122: Fluorapatite pattern in blue

### Step 3: Performing the Tune Cell Operation

1. Select the Fluorapatite pattern in the data tree.
2. Click **Tune Cell** in the Data Command panel  
— or —  
Right-click the pattern, then click **Tune Cell** on the context menu.  
The Tune cell dialog box will be displayed.
3. In the list of parameters, click **a**, and then click **h00** in the **Show** drop-down list to see lines which position depend only on **a**.



4. Zoom in on one of the (h00)-peaks.
5. The ghosts of pattern sticks are displayed in the graphical view. Use the slider to fine-tune **a** in order to fit the peak.




Fig. 123: (300) stick (located at  $2\theta=33.129^\circ$ ) of the pattern PDF 00-015-0876 shifted to fit the corresponding peak.

6. Check whether all peaks fit the pattern.
7. Repeat steps 3 through 6 for **c** by selecting **00l** in the **Show** drop-down list.
8. Use the **Replace** or **Append** buttons to either replace the current pattern or append the modified one to the data tree.

### Step 4: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name. (To overwrite an .EVA file, use the same name).
4. Click **Save**.




Fig. 124: Saving Francolite.EVA document

## Displaying a Pattern with Another Wavelength

The following procedure describes how to display a pattern with another wavelength.

The document used is held as a tutorial file, Francolite.RAW, found in the Tutorial directory.

The measured sample is a Carbonate-fluoroapatite (its mineral name is Francolite). We are going to check the pollution of the X-ray tube is suspected.

Steps:

---

1. Creating a new EVA document and importing Francolite.RAW.
  2. Performing the Search/Match operation.
  3. Changing the wavelength.
  4. Saving.
- 



### NOTE

A reference database is required to perform a Search/Match operation.

This tutorial was prepared using PDF 4+ 2009 as a reference database. You may have another reference database.

If you have no database, you cannot perform this part of the tutorial using the Search/Match.

---

## Step 1: Creating a New EVA Document and Importing Francolite.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed. Search the Tutorials/EVA<sup>1</sup> directory and select the Francolite.RAW file.
3. Click **Open**. The scan Francolite will be displayed in the graphical view of the EVA document.




Fig. 125: Import a Scan File dialog box




Fig. 126: Francolite.RAW imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

### Step 2: Performing the Search/Match Operation

Keep the default search parameters and run a search. Here are the results obtained:

Database PLU2009: 291440 - After Filters: 213504

Chemical Filter Database Filter Candidate List Selected Candidates

Index #	%	Source	ID	Q	S	I/Icor	Mineral	Name
<input checked="" type="checkbox"/>	1	86	PLU2009 PDF 00-031-0267	B	P		Yes	Carbonatefluorapatite
<input checked="" type="checkbox"/>	2	86	PLU2009 PDF 00-021-0141	B	D		Yes	Carbonate-fluorapatite
<input type="checkbox"/>	3	15	PLU2009 PDF 00-006-0454	B	D		Yes	Fluorellestadite, phosphatian
<input type="checkbox"/>	4	37	PLU2009 PDF 00-015-0876	S	P	1.5	Yes	Fluorapatite, syn
<input type="checkbox"/>	5	29	PLU2009 PDF 00-056-1014	O	P		No	Calcium Europium Fluoride Phosphate
<input type="checkbox"/>	6	60	PLU2009 PDF 00-047-0262	O	P		No	Sodium Calcium Hydrogen Carbonate Phosphate Hydrate

Group Duplicates Matched 114935 / 213504 Candidates in 15.3 s.

Carbonatefluoroapatite is easily identified. Select both patterns PDF 00-031-0267 and PDF 00-021-0141. It will allow viewing the change.




Fig. 127: Carbonatefluoroapatite pattern in blue. The two patterns are superposed.

### Step 3: Changing the Wavelength

1. Select the pattern PDF 00-021-0141 in the data tree.
2. In the pattern property table, clear the **Scan WL** check box to be able to enter another wavelength.

We want to check if there is a pollution of the X-ray tube. The  $WL\alpha_1$  radiation is the most frequently encountered spurious X-ray line in sealed X-ray tubes. Therefore, we are going to display the pattern for the  $WL\alpha_1$  radiation.

3. In the **Wavelength** drop-down list below, select the  $WL\alpha_1$  wavelength. The pattern PDF 00-021-0141 will be displayed with this new wavelength.
4. Zoom in from  $2\theta=20^\circ$  to  $2\theta=55^\circ$ .




Fig. 128: Carbonatefluorapatite pattern displayed with the  $WL\alpha_1$  wavelength in blue and original pattern in red. Zoom on the  $20-55^\circ$   $2\theta$  range.

5. The modified pattern looks to match some parts of the scan. It makes it possible to conclude there is likely a pollution of the X-ray tube.



### Step 4: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name. (To overwrite an .EVA file, use the same name).
4. Click **Save**.




Fig. 129: Saving Francolite\_WL.EVA document

## Working with the Waterfall View

The following procedure describes how to work with waterfall view.

The scan used is held as a tutorial file, Ep333f.RAW, found in the Tutorial directory.

Steps:

---

1. Creating a new EVA document and importing Ep333f.RAW.
  2. Creating and working with the waterfall view.
  3. Saving the EVA document containing the scan.
-

## Step 1: Creating a New EVA Document and Importing Ep333f.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed.
3. Search the Tutorials/EVA<sup>1</sup> directory and select the Ep333f.RAW file.
4. Click **Open**. The scan Ep333f will be displayed in the graphical view of the EVA document.




Fig. 130: Import a Scan File dialog box




Fig. 131: Ep333f.RAW file imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Creating and Working with the Waterfall View.

1. Select the **Waterfall display** check box in the Theta View property panel: the scans will be displayed as a waterfall.




Fig. 132: Waterfall display

2. Change the perspective of the display.

To change the perspective horizontally, click the X-axis to display the horizontal scroll bar. Point to the right end of the scroll bar. The mouse pointer will change into a double arrow. Move the double arrow left or right to reduce or expand the scroll bar.




Fig. 133: Perspective horizontally modified

To change the perspective vertically, proceed the same way but start by clicking on the Y-axis.




Fig. 134: Perspective horizontally and vertically modified

3. Zoom in on the right peak. The corresponding zone must be selected in the shadow representation at the back: move the cursor to one end of the targeted zoom area, press and hold the left mouse button while dragging the mouse until the cursor reaches the opposite end, then release the button.




Fig. 135: Zoom on the right peak

### Step 3: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name. (To overwrite an .EVA file, use the same name).
4. Click **Save**.




Fig. 136: Saving Ep333f.EVA document

## Working with the 2D View

The following procedure describes how to work with the 2D view.

The scan used is held as a tutorial file, Ep333f.RAW, found in the Tutorial directory.

Steps:

---

1. Creating a new EVA document and importing Ep333f.RAW.
  2. Creating and working with the 2D view.
  3. Saving the EVA document containing the scan.
-

## Step 1: Creating a New EVA Document and Importing Ep333f.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed.
3. Search the Tutorials/EVA<sup>1</sup> directory and select the Ep333f.RAW file.
4. Click **Open**. The scan Ep333f will be displayed in the graphical view of the EVA document.




Fig. 137: Import a Scan File dialog box




Fig. 138: Ep333f.RAW file imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.



## Step 2: Creating and working with the 2D view.

### Case #1: 2D view with intensity map

1. Make certain the list of scans is selected; if not, select it in the data tree.
2. Click **2D View** in the Create list of the Data Command panel  
— or —  
Right-click the scan list in the data tree to display the context menu. Click **Create** and then **2D View** on the related submenu.
3. The 2D View is displayed in a new tab. The Intensity map is displayed by default. Leave the **Intensity map** check box selected.




Fig. 139: 2D view with intensity map

4. Customize the intensity map by right-clicking the color scale on the right and choosing among the offered choices on the context menu. Select for example the **Diffraction Space Viewer** option.




Fig. 140: 2D view with colored intensity map

**Case #2: 2D view with levels**

1. Make certain the list of scans is selected; if not, select it in the data tree.
2. Click **Create Level** in the Tool list of the Data Command panel  
— or —

Right-click the scan list in the data tree to display the context menu. Click **Tool** and then **Create Level** on the related submenu.

The Create Level Theta dialog box will be displayed. By default 5 automatic levels will be created. They are displayed as ghost levels on the graphical view.




Fig. 141: Creating levels: Create Level dialog box and ghost levels in the graphical view

3. Click the **Append Levels** button and close the dialog box: the levels will be added in the data tree and displayed in the graphical view.




Fig. 142: Creating levels: levels displayed in the graphical view

4. Click **2D View** in the Create list of the Data Command panel  
— or —

Right-click the scan list in the data tree to display the context menu. Click **Create** and then **2D View** on the related submenu.

5. The 2D View will be displayed in a new tab. The Intensity map and the Levels are displayed by default. Leave the **Intensity map** and the **Levels** check boxes selected.




Fig. 143: 2D view with intensity map and levels

### Step 3: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name. (To overwrite an .EVA file, use the same name).
4. Click **Save**.




Fig. 144: Saving Ep333f-bis.EVA document

## Sorting a Multi-range Scan by Temperature in a 2D View

The following procedure describes how to sort scans by temperature from a non-ambient measurement in a 2D view.

The scan used is held as a tutorial file, Gui14.RAW, found in the Tutorial directory.

Steps:

---

1. Creating a new EVA document and importing Gui14.RAW.
  2. Creating the 2D view with levels.
  3. Sorting the scans on the Y-axis by temperature.
  4. Saving the EVA document containing the scan.
-

## Step 1: Creating a New EVA Document and Importing Guil4.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed.
3. Search the Tutorials/EVA<sup>1</sup> directory and select the Guil4.RAW file.
4. Click **Open**. The scan Guil4 will be displayed in the graphical view of the EVA document.




Fig. 145: Import a Scan File dialog box





Fig. 146: Guil4.RAW file imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Creating the 2D View with Levels

1. Make certain the list of scans is selected. If not, select it in the data tree.
2. Click **2D View** in the Create list of the Data Command panel  
— or —  
Right-click the scan list in the data tree to display the context menu. Click **Create** and then **2D View** on the related submenu.
3. The 2D View will be displayed in a new tab.



4. Click **Create Level** in the Tool list of the Data Command panel  
— or —  
Right-click the scan list in the data tree to display the context menu. Click **Tool** and then **Create Level** on the related submenu.  
The Create Level Theta dialog box will be displayed. By default 5 automatic levels will be created. They are displayed as ghost levels on the graphical view.




Fig. 147: Creating levels: Create Level dialog box

- Click the **Append Levels** button and close the dialog box: the levels will be added in the data tree and displayed in the graphical view.
- The Intensity map and the Levels are displayed by default. Leave the **Intensity map** and the **Levels** check boxes selected.




Fig. 148: 2D view with intensity map and levels, scans sorted by index on the Y-axis

### Step 3: Sorting the Scans on the Y-Axis by Temperature

- In the 2D View Property table, select **{SORT}** for the Left (Axis) property if necessary.
- In the 2D View Property table, select **Temperature** in the Y-axis drop-down list.




Fig. 149: 2D view with intensity map and levels, scans sorted by temperature on the Y-axis

### Step 4: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name. (To overwrite an .EVA file, use the same name).
4. Click **Save**.




Fig. 150: Saving Gui4.EVA document



## Working with the PIP and VIP Views

The following procedure describes how to create and manage PIP and VIP views.

The scan used is held as a tutorial file, m1.RAW, found in the Tutorial directory.

Steps:

---

1. Creating a new EVA document and importing m1.RAW.
  2. Creating and managing the PIP and VIP views.
  3. Saving the EVA document containing the scan.
-

## Step 1: Creating a New EVA Document and Importing m1.RAW



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed.
3. Search the Tutorials/EVA<sup>1</sup> directory and select the m1.RAW file.
4. Click **Open**. The scan m1 will be displayed in the graphical view of the EVA document.




Fig. 151: Import a Scan File dialog box




Fig. 152: m1.RAW file imported in the graphical view

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Creating and working with the PIP and VIP views.

### Case #1: Creating a PIP view



PIP mode button

1. Click the **PIP mode** button on the view toolbar  
— or —

Right-click anywhere in the graphical view to display the context menu. In the context menu, click the **PIP Mode** command.

A **PIP** text box will be added below the pointer.

2. Select the zone around the peak between the  $2\theta$  angles  $22^\circ$  and  $24^\circ$ : a window corresponding to the selected zone is created.




Fig. 153: Creating a PIP view

3. Move and resize the window to have a better view of the peak.




Fig. 154: Moving and resizing the PIP view

- Once you are satisfied with the PIP view, click anywhere in the graphical view. The PIP view is inserted in the graphical view and is linked to the corresponding zone in the scan.




Fig. 155: PIP view inserted and linked to the original zone

**Case #2: Creating a VIP view**

- Zoom in on the zone between the  $2\theta$  angles  $37^\circ$  and  $47^\circ$ .



VIP mode button

- Click the **VIP mode** button on the view toolbar  
— or —

Right-click anywhere in the graphical view to display the context menu. In the context menu, click the **VIP Mode** command. A **VIP** text box will be added below the pointer.

- Select the zone between the  $2\theta$  angles  $40.5^\circ$  and  $41.5^\circ$ : the VIP view editor is displayed.




Fig. 156: Creating a VIP view

4. Move the scale bar to increase the scale factor until 4.00 for example.




Fig. 157: Changing the scale

5. To insert the VIP view in the document, click anywhere in the graphical view.




Fig. 158: VIP view inserted in the graphical view

6. The VIP view can be modified by right-clicking anywhere in the graphical view, clicking **VIP mode** on the context menu and then **Edit**.

### Step 3: Saving

1. Click **Save As** on the **File** menu. The Save As dialog box will be displayed.
2. Select the appropriate disc drive and directory.
3. Type in the file name. (To overwrite an .EVA file, use the same name).
4. Click **Save**.




Fig. 159: Saving m1.EVA document

## Creating a Label from a Peak

Peaks can be used to create labels.

To do so:

1. Right-click a scan at the position you want to input a label (peak) to display the context menu.
2. Click **Create Peak at 2Th=...**
3. The peak is displayed in the graphical view and added to the peak list in the data tree.




Fig. 160: Peak inserted manually in the graphical view

4. Select the peak in the data tree if necessary.
5. In the Peak Property table, select the **Anchor lock** check box to be able to move the anchor and text position freely. Press the **Control** key and point to the scan to change the pointer into a hand. Move the hand to drag the anchor and the text. Clear the **Anchor lock** check box and proceed the same way to move the text position only.




Fig. 161: Anchor and text position changed

6. Clear the **Caption** text field and enter the desired text.




Fig. 162: Anchor and text position changed

7. Customize the resulting “label” if necessary:

- The peak color can be changed: in the Peak property table, select the desired color in the **Color** drop-down list. (Blue in the example).
- The background and text color can be changed: in the Peak property table, select the desired color in the **Background Color and Text Color** drop-down lists (Light gray and Dark violet in the example).
- The font size can be changed: in the Peak property table, enter the desired font size (in points) in the **Font Size** field (10 points in the example).
- The text rotation can be changed: in the Peak property table, enter desired the angle value in degrees in the **Text rotation** field ( $0^\circ$  in the example).
- The text margin can be changed: in the Peak property table, enter the desired text margin (in points) in the **Text margin** field (2 points in the example).
- The anchor can be removed or its style changed: in the Peak property table, select None or another anchor style in the **Anchor Style** drop-down list (Rectangle in the example).






Fig. 163: Label customized according to the procedure previously described

## Loading and Integrating Mergeable 2D Frames

The following procedure describes how to create an integrated scan from a set of frames which have been measured with different 2-theta angles to form a diagram like a Debye-Scherrer image.

The measurement used is stored as a tutorial file, "Cor standard.brml", found in the "Tutorial\2D Frame Data" directory.

### Steps:

---

1. Creating a new EVA document and importing Cor standard.brml.
  2. Selecting the integration cursor and drawing the integration area.
  3. Displaying the Cursors Preview tool and integrate.
-

## Step 1: Creating a New EVA Document and Importing the 2D Frames



New button



Import from file button

1. Click **New** on the **File** menu or use the toolbar button. The EVA document will be empty.
2. Click the **Import From File** data command or the **Import From File** button. The Import From File dialog box will be displayed.
3. Search the "Tutorials\EVA"<sup>1</sup> directory and select the "Cor standard.brml" file in the "2D Frame Data" sub-directory.
4. Click **Open**. The frames which were measured will be displayed in the frame view of the EVA document. A **Mergeable Frame List** will appear in the tree.




Fig. 164: 2D Frame View after loading a measurement with three mergeable frames

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

## Step 2: Selecting the Integration Cursor and Drawing the Integration Area



**Slice Cursor**  
button

1. Click on the **Slice Cursor** button.
2. Click with the left mouse button in the top left of the area which is to be integrated.
3. Drag the mouse to the bottom right until the desired integration area has been displayed. Release the left mouse button. The cursor will be added as child to the frame list node in the tree.
4. If necessary, adjust the integration area by clicking and dragging the handles of the cursor.




Fig. 165: Slice Cursor prepared for the subsequent integration

### Step 3: Displaying the Cursors Preview Tool and Integration

1. Click on the **Cursors Preview** command in the command bar or the context menu of the slice cursor.
2. If the preview is satisfactory, click on the Integrate button to create the final integrated scan. The scan is created in a separate scan list.




Fig. 166: The Cursors Preview Tool with the preview of the integration result




Fig. 167: The resulting integrated scan

## Loading and Integrating Stackable 2D Frames

The following procedure describes how to create an integrated scan from a set of frames which were measured with the same 2-theta angle and different theta angles to form a stackable frame list.

The measurement used is stored as a set of tutorial files, "HTSR258\*.gfrm", found in the "Tutorial\2D Frame Data\HTSR258" directory.

### Steps:

---

1. Creating a new EVA document and importing the 2D frames.
  2. Selecting the integration cursor and drawing the integration area.
  3. Displaying the Cursors Preview tool and integrate.
-

## Step 1: Creating a New EVA Document and Importing the 2D Frames



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document will be empty.
2. Click the **Import From File** data command or the **Import From File** button. The Import From File dialog box will be displayed.
3. Search the “Tutorials\EVA\2D Frame Data” directory and select all of the “HTSR258\*.gfm” files in the “HTSR258” sub-directory.
4. Click **Open**. The frames will be displayed in the frame view of the EVA document. A **Stackable Frame List** will appear in the tree.




Fig. 168: 2D Frame View after loading a measurement with stackable frames

## Step 2: Selecting the Integration Cursor and Drawing the Integration Area



**Ring Cursor**  
button

1. Click on the **Ring Cursor** button.
2. Click with the left mouse button in the left of the area which is to be integrated.
3. Drag the mouse to the right until the desired integration area has been displayed. Release the left mouse button. The cursor will be added as child to the frame list node in the tree.
4. If necessary, adjust the integration area by clicking and dragging the handles of the cursor.




Fig. 169: Ring Cursor prepared for the subsequent integration



### Step 3: Displaying the Cursors Preview Tool and Integrate

1. Change the direction of the integration in the cursor's properties to "2-Theta integration". Select "Merge as Single Scan".
2. Click on the **Cursors Preview** command in the command bar or the context menu of the ring cursor.
3. If the preview is satisfactory, click on the Integrate button to create the final integrated scan. The scan will be created in a separate scan list.




Fig. 170: The integration parameters of the ring cursor




Fig. 171: The Cursors Preview Tool with the preview of the integration result




Fig. 172: The resulting integrated scan

## Rocking Curve Integration on Stackable 2D Frames

The following procedure describes how to create an integrated scan from a set of frames which have been measured with the same 2-theta angle and different theta angles to form a stackable frame list.

The measurement used is stored as a set of tutorial files, "GM\_XRR.gfrm", found in the "Tutorial\2D Frame Data\GM\_XRR.gfrm" directory.

### Steps:

---

1. Creating a new EVA document and importing the 2D frames.
  2. Selecting the integration cursor and drawing the integration area.
  3. Displaying the Cursors Preview tool and integrate.
-

## Step 1: Creating a New EVA Document and Importing the 2D Frames



New button



Import from file button

1. Click **New** on the **File** menu or use the toolbar button. The EVA document will be empty.
2. Click the **Import From File** data command or the **Import From File** button. The Import From File dialog box will be displayed.
3. Search the “Tutorials\EVA\2D Frame Data” directory and select all of the “GM\_XRR\*.gfm” files in the “GM\_XRR” sub-directory.
4. Click **Open**. The frames will be displayed in the frame view of the EVA document. A **Stackable Frame List** will appear in the tree.




Fig. 173: 2D Frame View after loading a measurement with stackable frames

## Step 2: Selecting the Integration Cursor and Drawing the Integration Area



Area Cursor  
button

1. Zoom into the frame view by clicking and dragging the mouse to make the bright spots clearly visible.
2. Click on the **Area Cursor** button.
3. Click with the left mouse button in the top left of the area which is to be integrated.
4. Drag the mouse to the bottom right until the desired integration area has been displayed. Release the left mouse button. The cursor will be added as child to the frame list node in the tree.
5. If necessary, adjust the integration area by clicking and dragging the handles of the cursor.




Fig. 174: Ring Cursor prepared for the subsequent integration

### Step 3: Displaying the Cursors Preview Tool and Integrate

1. Change the X-axis of the integration in the cursor's properties to "Theta".
2. Click on the **Cursors Preview** command in the command bar or the context menu of the slice cursor.
3. If the preview is satisfactory, click on the Integrate button to create the final integrated scan. The scan will be created in a separate scan list.



Fig. 175: The integration parameters of the ring cursor




Fig. 176: The Cursors Preview Tool with the preview of the integration result




Fig. 177: The resulting integrated scan

## Adjusting the Default Mask of a 2D Frame

The following procedure describes how to adjust the default mask on a set of frames which have been measured with different 2-theta angles to form a diagram like a Debye-Scherrer image. Excluding the non-exposed area will lead to a more precise integration result.

The measurement used is stored as a tutorial file, "Cor standard.brml", found in the "Tutorial\2D Frame Data" directory.

### Steps:

---

1. Creating a new EVA document and importing Cor standard.brml.
  2. Check the integration with the Full Frame cursor.
  3. Changing the default mask to exclude non-exposed areas.
  4. Selecting the integration cursor and integrate.
-

## Step 1: Creating a New EVA Document and Importing the 2D Frames



New button



Import from file button

1. Click **New** on the **File** menu or use the toolbar button. The EVA document will be empty.
2. Click the **Import From File** data command or the **Import From File** button. The Import From File dialog box will be displayed.
3. Search the "Tutorials\EVA" directory and select the "Cor standard.brml" file in the "2D Frame Data" sub-directory.
4. Click **Open**. The frames which have been measured will be displayed in the frame view of the EVA document. A **Mergeable Frame List** will appear in the tree. The default masks which are part of the frame data will be drawn in green.




Fig. 178: 2D Frame View exhibiting non-exposed areas around the frames

## Step 2: Check the Integration with the Full Frame Cursor



Full Frame  
Cursor button

1. Click on the **Full Frame Cursor** button in the command bar. The whole merged frame is marked for integration.
2. Click the **Integrate Cursor** command in the command bar or the context menu of the cursor to create the integrated scan. The scan will be created in a separate scan list. A problem is clearly visible at the beginning of the scan where the intensity drops.





Fig. 179: Integrated scan with intensity drop at low angles due to mask misalignment




### Step 3: Changing the Default Mask to Exclude Non-Exposed Areas


1. Click on the **Edit Multiple Properties** check mark to display the combined properties for all three frames.



2. Click the **Use Custom Mask** check mark to display the mask properties:



3. Change the Radius to 980 pixel:



The displayed masks change according to the chosen radius. No non-exposed areas remain visible.




Fig. 180: 2D Frame View with an adjusted default mask

## Step 4: Selecting the Integration Cursor and Integration



Full Frame  
Cursor button

1. Click on the **Full Frame Cursor** button in the command bar. The whole merged frame is marked for integration.
2. Click the **Integrate Cursor** command in the command bar or the context menu of the cursor to create the integrated scan. The scan is created in a separate scan list.




Fig. 181: 2D Frame View with the full frame cursor drawn along its edges




Fig. 182: The resulting integrated scan

## Using Printing Options

This chapter will describe how to use some printing options.

### Vertical Table

“Vertical Table” is an option which allows “permuting” columns and rows. It proves useful when printing a column view with few objects but numerous properties.



New button



Import from file button

1. Click **New** on the **File** menu or use the button in the toolbar. The EVA document will be empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed.
3. Search the Tutorials/EVA<sup>1</sup> directory and select the ET20.RAW, ET50.RAW and ET80.RAW files.
4. Click **Open**. The scans will be displayed in the graphical view of the EVA document and added to the data tree.




Fig. 183: Import a Scan File dialog box




Fig. 183: ET20.RAW, ET50.RAW and ET80.RAW added to the data tree

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.

5. Select the scan list (2Theta list in the example).
6. Click **Scan Column View** in the Create list of the Data Command panel  
— or —

Right-click the scan list in the data tree to display the context menu. Click **Create** and then **Scan Column View** on the related submenu.

A Scan Column View will be created.

Visible	Icon	Color	Index	Name	Parent	Sample Name	File Name	Scan Type	Scan Status	Start
<input checked="" type="checkbox"/>			1	ET20.raw #1	2Theta	20% Boehmite/80% Corundum	ET20.raw	Coupled TwoTheta/Theta	Completed	22.000
<input checked="" type="checkbox"/>			2	ET50.raw #2	2Theta	50% Boehmite/50% Corundum	ET50.raw	Coupled TwoTheta/Theta	Completed	22.000
<input checked="" type="checkbox"/>			3	ET80.raw #3	2Theta	80% Boehmite/20% Corundum	ET80.raw	Coupled TwoTheta/Theta	Completed	22.000

7. Configure the table columns if necessary. To do so, right-click a column header to display the context menu. You can also drag and drop the columns to change their order.
8. In the Scan Column view Property table, select the **Vertical Table** check box. No change is visible in the scan column view. The resulting table will be visible in the Print preview.
9. Check the view is printable.
10. Click **Print Preview** in the View Data Command panel,  
— or —  
Right-click the scan column view in the data tree to display the context menu. Click **Print preview**. The print preview with the vertical table will be displayed.

2Theta

Visible	Yes	Yes	Yes
Icon			
Color			
Index	1	2	3
Name	ET20.raw #1	ET50.raw #2	ET80.raw #3
Parent	2Theta	2Theta	2Theta
Sample Name	20% Boehmite/80% Corundum	50% Boehmite/50% Corundum	80% Boehmite/20% Corundum
File Name	ET20.raw	ET50.raw	ET80.raw
Scan Type	Coupled TwoTheta/Theta	Coupled TwoTheta/Theta	Coupled TwoTheta/Theta
Scan Status	Completed	Completed	Completed
Start	22.000	22.000	22.000
End	45.000	45.000	45.000
Step Size	0.020	0.020	0.020
Time per Step	2.0	2.0	2.0
Temperature	25 (Room)	25 (Room)	25 (Room)
Time Started	0.0	0.0	0.0
Goniometer radius	200.5	200.5	200.5
2-theta	22.000	22.000	22.000
Theta	11.000	11.000	11.000
Chi	0.00	0.00	0.00
Phi	0.00	0.00	0.00
X-Drive			
Y-Drive			
Z-Drive			
Aux1	0.0	0.0	0.0

Fig. 184: Vertical table visible in the print preview

## Repeat Columns

“Repeat Columns” is an option which allows repeating the N first columns of a column view on each “row”. It proves useful when printing a column view with numerous objects and properties.

The “Repeat Columns” feature is available up to version 3.1



New button



Import from file button

1. Click **New** on the **File** menu or use the dedicated button in the toolbar. The EVA document is empty.
2. Click the **Import from file** data command or the **Import from file** button. The Import from files dialog box will be displayed.
3. Search the Tutorials/EVA<sup>1</sup> directory and select the fifteen first .RAW files from Aa.raw to Corundum-3.Raw in this example.
4. Click **Open**. The scans will be displayed in the graphical view of the EVA document and added to the data tree.





Fig. 185: Import a Scan File dialog box

<sup>1</sup> The directory path to the EVA tutorial files is "C:\Documents and Settings\All Users\Application Data\Bruker AXS\Tutorials" under Windows XP and "C:\ProgramData\Bruker AXS\Tutorials" under Windows 7. To access the folder with the explorer, select "Display hidden files and folders" in the Folders options.



Data	Description
Document	
Views	
2Theta View	Corund (Coupled TwoTheta/Theta)
Settings	1 Chemical Filter - 1 Database Filter
2Theta	14 Scans
Aa.raw #1	A-XRD-QU40DOL10COR20CAL30 (Coupled TwoTheta/Theta)
Alreftt.raw #2	ALREFLT (Coupled TwoTheta/Theta)
Asnosoll.raw #3	ASPIRIN - NO PRIMARY SOLLER SLIT (Coupled TwoTheta/Theta)
Assoller.raw #4	ASPIRIN - PRIMARY SOLLER SLITS (Coupled TwoTheta/Theta)
AUST_BC.raw #5	SAMPLE BC (Coupled TwoTheta/Theta)
AUST_EB.raw #6	SAMPLE EB (Coupled TwoTheta/Theta)
AUST_H.raw #7	SAMPLE H (Coupled TwoTheta/Theta)
Bb.raw #8	B-XRD-QU25DOL20COR20CAL35 (Coupled TwoTheta/Theta)
Bchips.raw #9	BCHIPS (Coupled TwoTheta/Theta)
BX100.raw #10	BX100 (Coupled TwoTheta/Theta)
Cc.raw #11	C-XRD-QU10DOL30COR20CAL40 (Coupled TwoTheta/Theta)
Corundum-1.raw #12	Corund (Coupled TwoTheta/Theta)
Corundum-2.raw #13	Corund (Coupled TwoTheta/Theta)
Corundum-3.raw #14	Corund (Coupled TwoTheta/Theta)

Fig. 186: Scans added to the data tree

5. Select the scan list (2Theta list in the example).
6. Click **Scan Column View** in the Create list of the Data Command panel  
— or —  
Right-click the scan list in the data tree to display the context menu. Click **Create** and then **Scan Column View** on the related submenu.  
A Scan Column View will be created.

Visible	Icon	Color	Index	Name	Parent	Sample Name	File Name	Scan Type
<input checked="" type="checkbox"/>		Red	1	Aa.raw #1	2Theta	A-XRD-QU40DOL10COR20CAL30	Aa.raw	Coupled TwoTheta/Theta
<input checked="" type="checkbox"/>		Blue	2	Alreftt.raw #2	2Theta	ALREFLT	Alreftt.raw	Coupled TwoTheta/Theta
<input checked="" type="checkbox"/>		Lime	3	Asnosoll.raw #3	2Theta	ASPIRIN - NO PRIMARY SOLLER SLIT	Asnosoll.raw	Coupled TwoTheta/Theta
<input checked="" type="checkbox"/>		Magenta	4	Assoller.raw #4	2Theta	ASPIRIN - PRIMARY SOLLER SLITS	Assoller.raw	Coupled TwoTheta/Theta
<input checked="" type="checkbox"/>		Black	5	AUST_BC.raw #5	2Theta	SAMPLE BC	AUST_BC.raw	Coupled TwoTheta/Theta
<input checked="" type="checkbox"/>		DarkRed	6	AUST_EB.raw #6	2Theta	SAMPLE EB	AUST_EB.raw	Coupled TwoTheta/Theta
<input checked="" type="checkbox"/>		DarkOr...	7	AUST_H.raw #7	2Theta	SAMPLE H	AUST_H.raw	Coupled TwoTheta/Theta
<input checked="" type="checkbox"/>		Green	8	Bb.raw #8	2Theta	B-XRD-QU25DOL20COR20CAL35	Bb.raw	Coupled TwoTheta/Theta
<input checked="" type="checkbox"/>		Cyan	9	Bchips.raw #9	2Theta	BCHIPS	Bchips.raw	Coupled TwoTheta/Theta
<input checked="" type="checkbox"/>		Yellow	10	BX100.raw #10	2Theta	BX100	BX100.raw	Coupled TwoTheta/Theta
<input checked="" type="checkbox"/>		LightSt...	11	Cc.raw #11	2Theta	C-XRD-QU10DOL30COR20CAL40	Cc.raw	Coupled TwoTheta/Theta
<input checked="" type="checkbox"/>		LightPink	12	Corundum-1.raw #12	2Theta	Corund	Corundum-1.raw	Coupled TwoTheta/Theta
<input checked="" type="checkbox"/>		Red	13	Corundum-2.raw #13	2Theta	Corund	Corundum-2.raw	Coupled TwoTheta/Theta
<input checked="" type="checkbox"/>		Blue	14	Corundum-3.raw #14	2Theta	Corund	Corundum-3.raw	Coupled TwoTheta/Theta

Fig. 187: Scan column view

7. In this example, it is interesting to configure the table columns differently. To do so, right-click a column header to display the context menu and access the configuration tools. You can also drag and drop the columns to change their order.

Name	Color	Sample Name	File Name	Start	End	Step Size	Time per Step	Time Started	Goniometer radius
Aa.raw #1	Red	A-XRD-QU40DOL10COR20CAL30	Aa.raw	23.990	46.990	0.020	2.0	0.0	200.5
Alreftt.raw #2	Blue	ALREFLT	Alreftt.raw	5.000	70.000	0.025	2.0	0.0	170.0
Asnosoll.raw #3	Lime	ASPIRIN - NO PRIMARY SOLLER SLIT	Asnosoll.raw	5.000	35.000	0.020	2.0	0.0	200.5
Assoller.raw #4	Magenta	ASPIRIN - PRIMARY SOLLER SLITS	Assoller.raw	5.000	35.000	0.020	2.0	0.0	220.5
AUST_BC.raw #5	Black	SAMPLE BC	AUST_BC.raw	20.000	40.000	0.020	3.7	0.0	200.5
AUST_EB.raw #6	DarkRed	SAMPLE EB	AUST_EB.raw	20.000	40.000	0.020	3.7	0.0	200.5
AUST_H.raw #7	DarkOr...	SAMPLE H	AUST_H.raw	20.000	40.000	0.020	3.7	0.0	200.5
Bb.raw #8	Green	B-XRD-QU25DOL20COR20CAL35	Bb.raw	23.990	46.990	0.020	2.0	0.0	200.5
Bchips.raw #9	Cyan	BCHIPS	Bchips.raw	5.000	70.000	0.025	2.0	0.0	170.0
BX100.raw #10	Yellow	BX100	BX100.raw	4.986	79.986	0.030	3.0	0.0	200.5
Cc.raw #11	LightSt...	C-XRD-QU10DOL30COR20CAL40	Cc.raw	23.990	46.990	0.020	2.0	0.0	200.5
Corundum-1.raw #12	LightPink	Corund	Corundum-1.raw	20.000	45.776	0.020	9.9	0.0	0.0
Corundum-2.raw #13	Red	Corund	Corundum-2.raw	45.797	71.573	0.020	19.8	151.8	0.0
Corundum-3.raw #14	Blue	Corund	Corundum-3.raw	71.593	97.370	0.020	39.6	441.6	0.0

Fig. 188: Scan column view with a new column configuration: file name and color in first

8. This table has many columns. Therefore it will not be possible to have all columns on each row when printing. It may be nice to keep, for example, the file name and its corresponding color at the beginning of each row.
9. In the Scan Column view Property table, enter 2 for the **Repeat Columns** value. It means the two first columns will be repeated on each row when printing.  
No change is visible in the scan column view. The resulting table will be visible in the Print preview.
10. Check the view is printable.
11. Click **Print Preview** in the View Data Command panel,  
— or —  
Right-click the scan column view in the data tree to display the context menu. Click **Print preview**.

The print preview will be displayed.

Print Preview(A4)

File PDFCreator 100%

### 2Theta

Name	Color	Sample Name	File Name	Start	End
Aa.raw#1	Red	A-XRD-QU40DOL10COR20CAL30	Aa.raw	23.990	46.990
Alreflt.raw#2	Blue	ALREFLT	Alreflt.raw	5.000	70.000
Asnosoll.raw#3	Green	ASPIRIN - NO PRIMARY SOLLER SLIT	Asnosoll.raw	5.000	35.000
Assoller.raw#4	Purple	ASPIRIN - PRIMARY SOLLER SLITS	Assoller.raw	5.000	35.000
AUST_BC.raw#5	Black	SAMPLE BC	AUST_BC.raw	20.000	40.000
AUST_EB.raw#6	Brown	SAMPLE EB	AUST_EB.raw	20.000	40.000
AUST_H.raw#7	Orange	SAMPLE H	AUST_H.raw	20.000	40.000
Bb.raw#8	Dark Green	B-XRD-QU25DOL20COR20CAL35	Bb.raw	23.990	46.990
Bchips.raw#9	Cyan	BCHIPS	Bchips.raw	5.000	70.000
BX100.raw#10	Yellow	BX100	BX100.raw	4.986	79.986
Cc.raw#11	Light Blue	C-XRD-QU10DOL30COR20CAL40	Cc.raw	23.990	46.990
Corundum-1.raw#12	Pink	Corund	Corundum-1.raw	20.000	45.776
Corundum-2.raw#13	Red	Corund	Corundum-2.raw	45.797	71.573
Corundum-3.raw#14	Blue	Corund	Corundum-3.raw	71.593	97.370

Name	Color	Step Size	Time per Step	Time Started	Goniometer radius	2-theta
Aa.raw#1	Red	0.020	2.0	0.0	200.5	23.990
Alreflt.raw#2	Blue	0.025	2.0	0.0	170.0	5.000
Asnosoll.raw#3	Green	0.020	2.0	0.0	200.5	5.000
Assoller.raw#4	Purple	0.020	2.0	0.0	220.5	5.000
AUST_BC.raw#5	Black	0.020	3.7	0.0	200.5	20.000
AUST_EB.raw#6	Brown	0.020	3.7	0.0	200.5	20.000
AUST_H.raw#7	Orange	0.020	3.7	0.0	200.5	20.000
Bb.raw#8	Dark Green	0.020	2.0	0.0	200.5	23.990
Bchips.raw#9	Cyan	0.025	2.0	0.0	170.0	5.000
BX100.raw#10	Yellow	0.030	3.0	0.0	200.5	4.986
Cc.raw#11	Light Blue	0.020	2.0	0.0	200.5	23.990
Corundum-1.raw#12	Pink	0.020	9.9	0.0	0.0	20.000
Corundum-2.raw#13	Red	0.020	19.8	151.8	0.0	45.797
Corundum-3.raw#14	Blue	0.020	39.6	441.6	0.0	71.593

Page 1 of 6 100%

Fig. 189: The two first columns "Name" and "Color" are repeated in the print preview