



NMR Training


- 1) How and where to sign-up for instrumentation time <http://argenta2.rd.unr.edu>
 - a. go over time rules listed on argenta2
 - b. make sure netid is prefixed with `unr\` when logging in to the instrument reservation system
- 2) Go over what nuclei can typically be ran on each NMR
 - a. 400-MR $^1\text{H}/^{19}\text{F}/^{31}\text{P}$
 - b. V400/V500 $^1\text{H}/^{19}\text{F}/^{13}\text{C}$
- 3) Go over magnet safety
 - a. No magnetic items in the NMR room
 - b. No ferromagnetic in NMR room
 - c. Keep items like wallets and cellphone at least six feet away from the magnet
 - d. No iPhone/ipads in the room during Helium fills (it has been shown that helium can penetrate the system clocks and render the devices dead)
 - e. Inform the NMR staff if you have medical implants
- 4) Putting a sample in a turbine and correctly using the depth gauge
 - a. Make sure you have at least three fingers of solvent or 760 microliters
- 5) Putting the sample/turbine into the magnet
 - a. switch on top of magnet
 - b. eject or insert button in VNMRJ parameters panel
 - c. e or i in the command line
- 6) Go over the parts of the VNMRJ interface
 - a. Experiments you'll typically use are: Proton, Carbon, gCOSY, ROESY, TOCSY, gHSQCAD, gHMBCAD
- 7) How to clear the automation queue in VNMRJ – makes a new folder to store data in
 - a. Automation -> Background Acquisition -> New Background Run
- 8) How to load in the default set of shims in VNMRJ *do this every time you run the NMR*
 - a. Tools -> Locator... ->double click on walkup_lk_d2o
 - b. Close locator window and then type su followed by enter into the command line
- 9) How to tune the probe on the 400-mr and v500: please note never to force the knob or barrel of the tuning stick. Tuning is a process of setting the NMR probe to the same frequency and impedance of the spectrometer
 - a. Tune the barrel on the capacitor stick to adjust the frequency, turn the knob to adjust the impedance; note the barrel and knob are interactive
 - b. In the command line type mtune
 - c. Select the rf channel and nucleus you want to tune
 - i. Channel 1: ^1H or ^{19}F
 - ii. Channel 2: ^{13}C , ^{31}P , ^{15}N
 - d. Click Start Probe tune button
 - e. Click Autoscale button
 - f. Adjust tuning knob and barrel on probe to get dip at correct frequency
 - g. Click Stop Probe tune button
 - h. Click Quit button
- 10) How to setup a ^1H 1D proton experiment: including when and why to change nt and d1
 - a. New study then fill out sample name (letters, numbers, underscore, and dash only) and select the deuterated solvent you're working in
 - b. Double click on the experiment you want, for example Proton, in the Study Queue and then adjust nt and d1
 - c. nt is adjusted based off concentration (about 1 to 5 mM for ^1H , 50 to 100 mM for ^{13}C for quick runs)
 - d. d1 is adjusted if your integrations do not come out to be whole numbers
- 11) How to submit a sample to the background
- 12) How to run multiple experiments on one sample
- 13) How to delete an experiment before submission




- 14) How to assess the quality of the shims (nl, dres on TMS or residual solvent proton)
 - a. the TMS peak should be a sharp symmetrical peak with a linewidth of less than 1 Hz nonspinning and less than 0.5 Hz spinning at 20 Hz.
- 15) Manual locking of the NMR
 - a. make sure Lock checkbox is unchecked
 - b. click lock scan
 - c. adjust z0 to be on resonance (sine wave becomes plateau)
 - d. adjust lock gain and phase so that the lock level value is stable and between 20 and 70 – a minor fluctuation of a few units is fine
 - e. do not adjust the lock phase
- 16) Manual shimming
 - a. generally, just adjust z1 through z7
 - b. if any shim is changes the value of the lock level significantly then start over from z1 again
- 17) How and where to save data to argenta
 - a. /mnt/argenta/mr400/data/PI/you
- 18) How to load data in Mestrenova

- a. click the open document button 
- b. navigate to your data folder on argenta
- c. open the FID file in your data folder


- 19) How to phase spectra in Mestrenova – ensuring the baseline is symmetrical around your peak

- a. click on the dropdown arrow on the phase button 
- b. select manual correction
- c. left and right drag in the purple box until the spectrum is properly phased

- 20) How to reference spectra in Mestrenova: both off TMS and the residual solvent proton

- a. click on the dropdown arrow of the reference button and select *Reference* 
- b. TMS peak – mouse over the peak until the peak turns red and click.
 - i. Type 0 in the *New Shift* box
- c. Residual solvent peak

- 21) How to integrate the peaks in Mestrenova

- a. click on the dropdown arrow of the integrate button and select *Manual* 
- b. drag over the peaks individually that you want to integrate
- c. set a specific integral to a known value
 - i. mouse over the integral until it glows, then right click and select Edit Integral...
 - ii. type in the new value in the **Normalized** box

- 22) How to label peaks with their chemical shifts in Mestrenova

- a. click on the dropdown arrow of the Peak Picking button and select *Peak by Peak*



- b. Mouse over the peaks that you want to label until the peak glows and then click

- 23) How to View tables for NMR data

- a. View → Tables...

- 24) How to measure coupling Constants

- a. Use crosshair tool and drag from the apex of one peak to the apex of another



- 25) How to predict chemical shifts based off a chemdraw structure in Mestrenova
 - a. Draw molecule in Chemdraw
 - b. Copy and paste molecule into a blank page in Mestrenova
 - c. In Mestrenova Predict → 1H Spectrum
- 26) How to change display parameters in Mestrenova: for example spectrum color, font size, etc.....
- 27) How to save processed spectra in Mestrenova
- 28) How to print processed spectra in Mestrenova
- 29) How to copy and paste spectra in Mestrenova
- 30) How to create your data folders on argenta
- 31) Fill out instrument problem sheet on argenta
- 32) No food or drink in the lab
- 33) Point out the oxygen sensor in the lab. If it goes off leave the room immediately
- 34) In the event the building alarm goes off. Leave the room immediately
- 35) Run three checkout samples: CHCl₃ in acetone (manual shim less than 0.5 Hz spinning, less than 1 Hz nonspinning), low volume D₂O, DMSO
- 36) How to manually run a sample
 - a. Forget about the study queue and don't use it
 - b. Load in the parameters for a Proton by clicking Proton in the available experiments
 - c. Tune the probe using the mtune macro
 - d. Manually lock the NMR
 - e. Manually shim the NMR
 - f. Adjust nt and d1
 - g. Type **ga** to go and acquire the data
 - h. Save your data to argenta
- 37) You must run your first five samples in the presence of the SIL staff
 - a. Sample 1 date: SS 1/30
 - b. Sample 2 date: SS 1/30
 - c. Sample 3 date: _____
 - d. Sample 4 date: _____
 - e. Sample 5 date: _____

Trainee Signature

Trainee Printed Name

Trainee Email: _____

Date Training completed: _____

Trainer Signature

Trainer Printed Name